**Operating Manual** 

# **Resting and Stress ECG**



Part 3: custo diagnostic software for custo cardio 300/400



MSW 0030 - DK 1834 Version 001 - 10/01/2020





© 2020 custo med GmbH

This Operating Manual may not be copied in its entirety or in part, duplicated in any form or by any means or translated into another language without the prior written consent of custo med GmbH.

The manufacturer reserves the right to change the information in this Operating Manual without prior notice. The current version can be downloaded from our website: www.customed.de.

#### CAUTION:

This Operating Manual is part of a modular system, consisting of four parts. All four parts must be downloaded from the Internet or from a CD to ensure the Operating Manual is complete.

## **Operating Manual**

## **Resting and Stress ECG**



Part 3: custo diagnostic software for custo cardio 300/400

## Table of contents

3.1	Symbo	Is used in this Operating Manual	4
3.2	custo c	liagnostic program structure	5
3.3	Connec	cting custo med devices to the PC	6
	3.3.1	Connecting and configuring custo cardio 300	6
	3.3.2	Connecting and configuring custo cardio 400 BT	8
	3.3.3	Connecting training devices for stress ECG	11
	3.3.4	Configuring a training device for stress ECG	12
3.4	ECG re	cording with custo diagnostic	13
	3.4.1	Recording of resting ECG	14
	3.4.2	Recording of stress ECG	17
	3.4.3	Extended ECG settings	24
3.5	Workin	ng with the evaluation	
	3.5.1	Opening an evaluation via the evaluation search	
	3.5.2	Opening an evaluation via the examination menu	
	3.5.3	Resting ECG evaluation	
	3.5.4	Stress ECG evaluation	
	3.5.5	Confirming the evaluation	
	3.5.6	Ending the evaluation	
	3.5.7	Optional: Reporting with approval process	
3.6	Append	dix	
	3.6.1	Values and formulas in the ECG evaluation	
	3.6.2	Keyboard navigation and shortcuts in custo diagnostic	41

## 3.1 Symbols used in this Operating Manual

ACTIONS THAT ARE PROHIBITED or not allowed under any circumstances!	$\bigcirc$
WARNING used to indicate situations which, if not avoided, could result in personal injury or property damage	
NOTE provides important information which must be observed	
TIP contains practical information to assist you with your work	
Words highlighted in colour indicate buttons or click paths to the corresponding program point, e.g. Examination, Settings,	Words highlighted in colour

### 3.2 custo diagnostic program structure

The custo diagnostic program is divided into three areas – User, Patient and Examination. This structure ensures that you can always recognise who (which user) is carrying out what type of examination with whom (which patient). The main menus of each area can be reached by clicking on User, Patient or Examination.

In the User **1** main menu, the user of the system can be selected. The administration of users takes place in the custo diagnostic service centre (creating users, allocating user rights, user-specific settings).

The **Patient 2** main menu is used for patient management. Its most important functions include Search for patient, Add new patient and Search Evaluation.

The Examination ③ main menu lists all of the examination types which are possible with custo diagnostic. Modules already purchased are active (black font), all other modules are inactive (light grey font).

This menu is also linked to the Settings **4** area. This area is for making crossprogram, examination-related and user-specific settings.



### 3.3 Connecting custo med devices to the PC

**Prerequisite:** custo diagnostic is installed on your PC and ready to use. The custo med devices and components may only be connected to the PC after custo diagnostic has been installed. The required device drivers are installed on the PC via the custo diagnostic standard setup or by specific selection during the custo diagnostic setup.

#### 3.3.1 Connecting and configuring custo cardio 300

#### custo cardio 300 USB version

Connect the ECG device to the PC using the supplied USB cable. Power is supplied via the USB cable. This charges devices with an integrated battery. The Windows driver installation starts automatically. After the driver installation has been completed, the device is configured in custo diagnostic:

- Start custo diagnostic and select Examination, Resting ECG or Stress ECG, Settings, Device 1, ECG Device 2.
- The ECG device is shown in the "ECG devices" section: cardio 300 USB 3.
- If it does not appear, click Scan 4.
- Select the ECG device cardio 300 USB 3.
- Later on, the ECG recording can also be started by pressing the
   button on the device. Here you can define the recording type:
   Auto start (3) (automatic ECG 10 s) or Start (3) (manual recording, any duration).
- Save 7 your input and close the screen with End 8. The ECG device is ready for operation.
- > For stress ECG: Configure the training device, see Chapter 3.3.3 and 3.3.4.

ECG Device       Network         ECG Devices       Settings for ecg device pump         Cardio 300 USB (SN: UBE00300)       Suction power (06)         Cardio 400 BT (SN: UBE00300)       Detachment time (060 sec)         Cardio 400 USB (SN: 00 0001 00000)       Detach electrodes at end of test auton         QRS Trigger       Scan Wifi         Add Cardm       Custo guard ECG         Data rate       1 ms (2 ms (4 ms (5 ms	Device	Export	lenu/Functions		<ul> <li>Print</li> </ul>	•		ECG	ng ECG	ng ECG	ECG	•	Print	rint	Print	Print	Print	Print	nt										M	en	u/ru	ncti	ons		Exp	ort						] [	Der	vic	:e	J				Jiag	gno	stic					4	•	
ECG Devices       Settings for ecg device pump            Cardio 300 USB (SN: UBE00300)         Cardio 300 BT (SN: UBE00300)         Cardio 400 BT (SN: UB 0001)         Cardio 400 USB (SN: 00 0001 00000)         Ceardio 400 USB (SN: 00 0001 00000)         Cleaning time (0180 min)         Detach electrodes at end of test auton         Custo router         Port 4001         Custo guard ECG         Data rate         1 ms 2 ms 4 ms 8 ms         ECG Device Button Action         At device button pressed         Same and a second			etwork	evice 2	ECG D								ECG D	CG De	ECG [	ECG	ECG	CG I	G De	De	De	De	Dev	ev	lice	•	2		Ne	etv	ork																										4	•	
Cardio 300 USB (SN: UBE00300)       Suction power (06)         Cardio 300 BT (SN: UBE00300)       Detachment time (060 sec)         Cardio 400 USB (SN: 00 0001 00000)       Cardio 400 USB (SN: 00 0001 00000)         Cardio 400 USB (SN: 00 0001 00000)       Detachment time (060 sec)         Cleaning time (0180 min)       Detach electrodes at end of test auton         Custo guard ECG       ECG Device Button Action         Data rate       1 ms (2 ms) 4 ms (8 ms)	s for ecg device pump	Settin						vices	Devices	Devices	vices																							l			5	Set	ttin	ngs	s f	or	ecç	g c	dev	vice	e pu	mp											
Cardio 300 BT (SN: UBE00300) Cardio 400 BT (SN: UB 0001) Cardio 400 USB (SN: 00 0001 00000) QRS Trigger Scan 4 Scan Wifi Add Cardm Custo router Port 4001 Custo guard ECG ECG Device Button Action At device button pressed 3	1 power (06)	Suctio			1: UBE00300)	B (SN: UB	USB (	5 300 US	dio 300	rdio 300	o 300 USI	B (SN: UE	BE00300	0300)	200300	E0030	2003	00300	300)	00)	0)	0)	))	)									-					Suc	ctic	on	p	owe	er	(0		6)						-	•	Ь	ev	el 3	3		1
□ cardio 400 USB (SN: 00 0001 00000)       □ Detectimient time (0180 min)         □ QRS Trigger       Scan       4 Scan Wifi         △ QRS Trigger       Scan       4 Scan Wifi         △ QRS Trigger       Scan       4 Scan Wifi         △ Custo router       Port       4001         Custo guard ECG       ECG Device Button Action         □ Data rate       1 ms       2 ms       4 ms       8 ms	mont time (0 60 coc)	Dotac			UBE00300)	(SN: UBE	BT (S	0 300 BT	dio 300	irdio 300	0 300 BT	(SN: UBE	000300)	300)	0300)	00300	0300	300) 11)	)0)	)	)	)											Г	1			١.	Det	tac	-br		+	tin	-		0	60										20		
QRS Trigger       Scan       4       Scan Wifi       Add Cardm         custo router       Port       4001         Custo guard ECG       ECG Device Button Action         Data rate       1 ms       2 ms       4 ms       8 ms	a time (000 sec)	Cloan		000)	0.00001	(SN. 0D C	UICE /	400 US	dio 400	rdio 400	0 400 US	(SN. 00 (	0001)	01.00	0001	001)	001)	001 (	1 00	0.01	000	000	000	00	00	22													an	un.		tim	~	(0	. (1	1 0 0	00 i	n)	, 							H	30		 21
QRS Trigger       Scan Wifi       Add Cardm         custo router       Port       4001         Custo guard ECG       ECG Device Button Action         Data rate       1 ms       2 ms       4 ms       8 ms	ach electrodes at end o					, on 00			400	100 100	0 100 001	5 (54.00	,			0001	0001	001 (	1 00	00						.,													De	ota	act act	h o		tre	od	00	at c	nd	oft	oct	aut	om	ati	cal	llv				 "
custo guard ECG     ECG Device Button Action       Data rate     1 ms     2 ms     4 ms     8 ms			Add Cardm	Scan Wifi	can 4	Scan	er	Trigger	S Trigge	RS Trigge	Trigger	Scan	4	4	4	4	4	4	4					5	Sca	an	Wif	fi	+	,	Add	Car	dm																										
custo guard ECG         ECG Device Button Action           Data rate         1 ms         2 ms         4 ms         8 ms           At device button pressed         §			1002	101					, outer	, outer																																																	
Data rate     1 ms  2 ms  4 ms  8 ms  At device button pressed	vice Button Action	ECG D					ECG	ard ECG	guard E	o guard E	uard ECG																							1			E	ECO	GC	De	vic	e I	But	tto	on	Act	tion												
	ce button pressed	At de	) 8 ms	○ 4 ms	ns 🔿 2 ms	) 1 ms (	٠	e 🖲	rate	rate	te 🖲	) 1 ms	○ 2 ms	2 ms	) 2 m	) 2 m	2	2 m	ms	ns	าร	าร	s	5	0	) •	4 m	ns	С	) 8	3 ms						'	At (	de	vio	ce	bu	tto	on	pr	ess	ed					5	0	Au Sti	ito art	Sti	art		
																			_	_	_	_	_	_				_	_	_	_	_	_		_		_	_	_		_	_	_	_		_		_	_	_	_	_	_		_		_		





## custo cardio 300 Bluetooth & USB version<sup>1)</sup> with rechargeable batteries or batteries

#### Setting up the Bluetooth connection

- Plug the Bluetooth USB stick into the PC. The driver installation starts automatically.
- Check whether the Bluetooth driver was installed correctly:
- On your Windows desktop right-click on Workspace or Computer.
   Select Manage in the context menu.
- In the left side of the window, click Device Manager.
- In the right side of the window, expand the Bluetooth item. The options Broadcom BCM20702 Bluetooth 4.0 USB Device and Microsoft Bluetooth Enumerator should appear.

#### Connecting custo cardio 300 BT to the PC

- Switch the device on by pressing and holding the custo med logo ( until the battery symbol, Bluetooth symbol and electrode monitoring field illuminate.
- > Open Windows Control Panel.
- > Click on Devices and printers, Add device. The ECG device is found.
- Select the custo cardio 300 entry in the "Add device" dialogue box and click on Next. The device is added without pairing code.

Device configuration for resting and stress ECG

- Start custo diagnostic and select Examination, Resting ECG or Stress ECG, Settings, Device 1, ECG Device 2.
- The ECG device is shown in the "ECG Devices" section: cardio 300 BT 3.
- ▶ If it does not appear, click Scan ④.
- Select the ECG device cardio 300 BT 3.
- Later on, the ECG recording can also be started by pressing the button on the device instead of using the software interface. Define the recording type:
   Auto start (automatic ECG 10 s) or Start (manual recording, any duration).
- Save vour input and close the screen with End vous. The ECG device is ready for operation.
- > For stress ECG: Configure the training device, see Chapter 3.3.3 and 3.3.4

Resting ECG	•	Print	Menu/Fu	nctions	Export	Device	Diagnostic	4	•
		ECG Device 2	Network					•	•
ECG Devices					Setting	s for ecg device pum	>		
ardio 300 USB	(SN: UB	E00300)		*	Suction	n power (06)	▼ Le	evel 3	
cardio 300 BT (9)	SN: UBE	00300) 1001)			Detach	ment time (0 60 sec	-)	30	sec
cardio 400 USB	(SN: 00	0001 00000)			Cleanir	na time (0180 min)	·/	30	mir
	(0.11.00	,			Det	tach electrodes at end	of test automatical	lv lv	
QRS Trigger	Scan	4 Scan Wifi	Add	Cardm					
custo router		F	Port 4001						
custo guard ECG					ECG D	evice Button Action			
Data rate 🛛 🔘	1 ms (	⊃2 ms ○4 m:	5 () 8 ms		At dev	ice button pressed	5 @ Au 6 _ Sta	to Start Irt	

🖸 custo cardio 300 BT&USB



1) To use the USB port, follow the instructions on the previous page.

## 3.3.2 Connecting and configuring custo cardio 400 BT (with Bluetooth and USB)

Connect the supplied power supply unit to the custo cardio 400 BT ECG device or to the custo cardio 400 BT supply line of the telescopic boom. Connect the power supply unit to the mains.

Either both connections (BT and USB) or only one connection can be configured as necessary. The user can define the configuration sequence.

#### Setting up the Bluetooth connection

- Plug the Bluetooth USB stick into the PC. The driver installation starts automatically.
- > Check whether the Bluetooth driver was installed correctly:
- On your Windows desktop right-click on Workspace or Computer.
   Select Manage in the context menu.
- In the left side of the window, click Device Manager.
- In the right side of the window, expand the Bluetooth item. The options Broadcom BCM20702 Bluetooth 4.0 USB Device and Microsoft Bluetooth Enumerator should appear.

Connecting custo cardio 400 BT Bluetooth to the PC

- > Open Windows Control Panel.
- Click on Devices and printers, Add device. The ECG device is found.
- Select the custo cardio 400 BT entry in the "Add device" dialogue box and click on Next. The device is added without pairing code.

#### Setting up the custo cardio 400 BT USB connection

Connect the ECG device to the PC using the supplied USB cable. The Windows driver installation for the new hardware starts automatically. Select the "Find driver automatically" option. After the driver installation has been completed, the device is configured in custo diagnostic.

The device configuration for the resting and stress ECG is on the next page.



Device configuration for resting and stress ECG

- Start custo diagnostic and select Examination, Resting ECG or Stress ECG, Settings, Device 1, ECG device 2.
- The ECG device is shown in the "ECG devices" section: cardio 400 BT and custo cardio 400 USB 3.
- If it does not appear, click Scan 4.
- Select cardio 400 BT and/or custo cardio 400 USB 3.
- Later on, the ECG recording can also be started by pressing the

   D button on the device instead of using the software interface.
   Define the recording type: Auto start () (automatic ECG 10 s)
   or Start () (manual recording, any duration).
- Save ② your input and close the screen with End ③. The ECG device is ready for operation.
- > For stress ECG: Configure the training device, see Chapter 3.3.3 and 3.3.4.



		User			custo med	GMDH			
		Patient	n		Resting E	CG			
Resting ECG	-	Print	м	lenu/Functions	Export	Device 1	Diagnostic	•	•
		ECG Device	2 N	etwork				•	•
ECG Devices					Settir	ngs for ecg device pur	ıp		
🗌 cardio 300 USB	(SN: UE	3E00300)		*	Suction	on power (06)	•	Level 3	-
cardio 300 BT (S cardio 400 BT (S	SN: UBE	00300)			Deta	hment time (060 se	ic)	30	sec
cardio 400 USB	(SN: 00	0001 00000)			Clean	ing time (0180 min)	)	30	min
					I D	etach electrodes at en	d of test automa	tically	
				~					
QRS Trigger	Scan	4 Scan V	Vifi	Add Cardm					
custo router			Port	4001					
custo guard ECG					ECG I	Device Button Action			
Data rate 🔹 🖲	1 ms	○ 2 ms ○ 4	ms 🤇	) 8 ms	At de	vice button pressed	5 e 6 ·	Auto Start Start	
								End	0

#### Additional custo cardio 400 BT device settings

The settings for the device pump (suction power, etc.) can also be found under Examination, Resting ECG or Stress ECG, Settings, Device **1**, ECG Device **2**.

#### **3** Suction power

Level 3 is the default setting for resting ECG, and level 5 is the default setting for stress ECG. You can change the suction power at a later time on the device or use the automatic suction level control.

With automatic suction power control, the suction level is increased until the system is airtight and the electrodes rest tightly enough against the skin of the patient. If the suction level is set manually, automatic suction level control is deactivated.

#### Oetachment time

Blow the moisture out of the electrode cables after an examination is "Detach electrodes at end of test automatically" is selected.

#### **5** Cleaning time

Extended blowing-out of moisture out of the electrode cables, e.g. at the end of a working day.

#### **6** Detach electrodes at end of test automatically

Selecting this option automatically detaches the electrodes from the patient after an examination and the length of time set under "Blowing after test" runs as specified. Disable this option in the resting ECG settings if a stress ECG is to be performed immediately following a resting ECG.

- > Click on Save (bottom left) to apply your settings.
- > Close the screen with End (bottom right).
- > The device is ready for operation.

		User		custo med	GmbH			? _
		Patient						
		Examination		Resting EC	G			
Resting ECG	-	Print	Menu/Functions	Export	Device 1	Diagnostic		
		ECG Device 2	Network					
ECG Devices				Setting	gs for ecg device pur	np		
🗌 cardio 300 USE	8 (SN: UB	E00300)		3 Suction	n power (06)	•	Level 3	•
🗌 cardio 300 BT	(SN: UBE	00300)						
🔳 cardio 400 BT	(SN: UB (	0001)		4 Detach	nment time (060 s	ec)	30	sec
🔳 cardio 400 USE	8 (SN: 00	0001 00000)		5 Cleanir	ng time (0180 min	)	30	min
				🌀 🔳 Det	tach electrodes at er	nd of test automat	ically	
			-					
QRS Trigger	Scan	Scan Wifi	Add Cardm					



#### 3.3.3 Connecting training devices for stress ECG

Connect the training device to the PC using the supplied cable (serial interface <sup>1</sup>). Make a note of the number of the serial connection, see Windows device manager. The number of the serial connection will be needed later in custo diagnostic.

Some devices (e.g. ergometer ec5000 and treadmill er2100) can be connected to the PC using a network cable.

Start custo diagnostic and open the Examination, Settings, Interface, Devices 1 screen. In the left side of the window, select the device, e.g. Ergometers, No. 1 2. In the right side of the window, change any settings as required. If the device type is known, select the Device 3 option and in the "Device" dropdown list select the device type, e.g. ec5000 4. If the device type is not obviously apparent, select the Protocol 3 option and in the "Protocol" dropdown list select the connected device, e.g. customed/ergoline 6.

In the "Options" dropdown list, set the device options **2**. Details about the device connection are entered in the "Interface" area **3**.

Use the Test ③ button to check whether the connection between device and PC is working. If the connection is working, the test dialogue box shows "Status: O.K. started".

Save 10 your input and close the screen with End 10.

		User		custo med G	mbH		? _ ×	ned GmbH		
		Patient								
		Examination		Resting ECG			•	) ECG		
System	-	Database	Export	System	Interface	< →		I	nterface	
		Devices 1	i.s. medical band	custo connect	com RF	< >		ect c	om RF	
Device se	tting for the work	station								
Deviles						•				•
Device	Ergometers		•	Settings made b	/ <b>3</b> Device	5 Protocol	_	ade by	O Device	5 • Pro
	No. 1, ec500	0 (SIO 30) 🗾			•				0	
	-No. 2, unkno	wn		Devi	ce 4 ec5000		-	Protocol	o customed/ergo	line
	No. 3, unkno	wn								
	No. 4, Unkno	wn		Opti	ons SpO2 data		<u>^</u>	Options	SpO2 data	
	No. 5, Unkno	wn			BP data				BP data	
	No. 6, unkno	wn			7 electr. sadd	lle height adjustment			7 electr. sadd	le height adjus
	No. 2, unkno	wii			ECG transm	hission			ECG transm	ISSION
	-No. 9, unkno	wn			Simulation				Simulation	
	No. 10 unkn	000			disable wat	chdog	·		disable wate	chdog
	No. 11 unkn	own		Interface 🔘	Serial/USB () TCP/I	P		Seria	I/USB () TCP/I	P
	-No. 12. unkn	own		8				8		
	No. 13. unkn	own			COM COM 30 -	Baudrate 2400	-		COM COM 30 A	Baudrate
	No. 14, unkn	own					_			
	No. 15. unkn	own		Description	ec5000		_		customed/ergo	line
	No. 16, unkn	own				Tect				т
						itse 🔰				

1) The number of serial interfaces on the PC can be expanded using USB-to-serial converters or a PCI plug-in card with serial connections.

#### 3.3.4 Configuring a training device for stress ECG

- In custo diagnostic, select
   Examination, Stress ECG, Settings, Device ①, Training device ②.
- > Select the Ergometer 3 or Treadmill 4 that was previously set.
- Select the blood pressure module of the ergometer or the previously connected and set sphygmomanometer (5).
- Select the SPO2 module of the ergometer or the previously connected and set SPO2 meter 6.
- Click on Save I to apply your input.
- Click on End <sup>3</sup> to close the screen.
- > The training device is ready for operation.

		User		custo med	GmbH		? _
		Patient					
		Examination		Stress ECG			
Stress ECG	•	Print	Menu/Functions	Export	Device 1	Diagnostic	< →
		ECG Device	Training Device	Network			< →
			2				
Ergometer							
-							
Ergometer •	NO. 1,	ec5000 (SIO -	Test	-			
Blood Pressure •	Ergom	neter ·	BP Control	-			
SPO2 -	Manua	al 🔺	Test	1			
Treadmill							
Treadmill	Manua	al 🔺	Tost				
Blood Pressure V	Manua	al 🔺	1600				
SPO2	Manus	al 🔺	Tost				
5102	Piuria		1650	_			
Treadmill accelerat	tion						
<ul> <li>Stress Level 3</li> </ul>		•	0.185 m/s <sup>2</sup>				
Save							enu 🧿

### 3.4 ECG recording with custo diagnostic

#### Note on the procedure

The steps necessary to record and evaluate ECG data in custo diagnostic are shown without a surgery IT system or HIS connection.

#### Program start, calling the ECG

- Start custo diagnostic and log in.
- Click on Examination, Resting ECG ①, New Resting ECG ② or: Examination, Stress ECG ③, New Stress ECG.

#### Selecting the patient

- Select a patient for the examination:
   Enter the patient's name into the input fields in the search screen ④.
- Select the patient from the list **5**.
   Confirm your selection with <u>Select Patient</u> **6**.
   You can also select the patient by double-clicking on the name.

#### New patient

- If the patient does not yet exist in your database:
   Click on New Patient 2.
- Enter the patient data.
   The fields marked with an asterisk are mandatory.
- > Save the data, the patient is entered into the database.

#### Further work steps

- > 3.4.1 Recording of resting ECG.....14

	User		User			User			custo med GmbH
	Patient		Patient			Patient			
	Examin		Examin			Examination			Resting ECG
Holter		New Resting ECG 2		Last name		Mustermann 4	)		
ABPM				First name					
Resting ECG 1		Show Evaluation		Patient ID					
Stress ECG 3		Show Comparison		Patient Group		All patients		•	
Cardiopulmonary Exercise Testing				Assignment	Physician	All physicians		-	
Spirometry		Assign Evaluation			Physician ID				
Telemetry									
Cardiac Rehab		Settings		Last name	First name	Date of birth	Pat. ID		
Mobile Cardiac Telemetry				Mustermann	Absoluta	10.10.1960	000000001	^	
Blood glucose				Mustermann	Franz 5	10.10.1960	0000000002		
Ecg Streaming									
Task Manager									
Worklist									. custo
								•	
							2 patient	(s)	
Settings				Select Patient	6	Edit Patient			
				New Patient	7				
Cancel		Cancel		Cancel					



#### 3.4.1 Recording of resting ECG

If several ECG devices are connected to the workplace, the "Select ECG Device" **1** dialogue will be displayed. Select the ECG device **2** and click on Confirm. If only one ECG device is connected, this step can be omitted.

Select ECG Device	2-593
cardio 300 USB (SN: UBE00300) cardio 300 BT (SN: UBE00300)	<u> </u>
Add Confirm 3 Cance	iagnostic

#### Monitoring and electrode control

The patient's ECG signal will be shown on the display but not recorded (monitoring). Work steps before the start:

- Change the type of lead if necessary 4.
- Check if all electrodes are attached optimally.
   If there are red lines on the screen, the contact between the skin and electrode(s) is insufficient. The corresponding electrodes will need to be reattached. For ECG devices with USB connection, selecting Electrodes (9) provides you with a graphic representation of the signal quality.
- > Set the required filter (Options menu 6).



The preset standard procedures for automatic ECG and manual recordings are described here. These procedures can be changed in the custo diagnostic settings, see Examination, Resting ECG, Settings, Menu/Functions, Workflow.

#### Note on recording with a tablet PC:

The rechargeable battery capacity of the system is queried before a recording is started. If it is less than 15%, no new recording can be performed.

#### Automatic ECG<sup>1)</sup> – Auto Start button 🥑

- Click on Auto Start to start the automatic recording. The default setting for the duration of automatic recording is ten seconds.
- At the end of the ten seconds, the recording is automatically ended, saved, measured and printed out.

#### Manual recording<sup>1)</sup> – Start <sup>(3)</sup> button or Enter key:

- If you would like to perform a recording without a time limitation, (e.g. because of suspected irregularities), trigger the recording using the Start ③ key.
- At least ten seconds of the ECG must be recorded before a recording can be ended.
- > Use Stop to end the recording, the ECG interface remains open.
- Use Start 8 and Stop to record additional sections.
- > Use End 9 to close the recording.
- Click on Confirm 1 in the End dialogue to ensure that the recording is saved, measured and displayed as an evaluation.







1) The preset standard procedures for automatic ECG and manual recordings can be changed in custo diagnostic under Examination, Resting ECG, Settings, Menu/Functions, Workflow. Save your input.



Editing options during the recording process

#### > Automatic marking of the ECG

Clicking on the Marker 1 button will automatically mark the last six seconds of the recording. A dialogue appears for specifying <sup>2</sup>), printing and saving the marked part 12.

#### Manual marking of the ECG

By clicking on the Start marking <sup>(1)</sup> button you can determine the length of the marking yourself. The marking continues to run until you click on the End marking button. A dialogue appears for specifying<sup>2)</sup>, printing and saving the marked part <sup>(2)</sup>.





Tip: Keyboard shortcut



#### 2) Text modules for specifying marked parts

Use: If a part is marked during the recording, the "Markings" dialogue appears. The marked parts will be automatically specified by pressing the corresponding key (e.g. F5) or a previously configured text module button.

Configuring text modules: Under Examination, Resting ECG, Setings,

Menu/Functions, Markings, the text modules can be configured for specifying marked parts. Enter a name for the text module in the "Description" field. This name later appears on the button for calling the text module. In the "Text module' field, enter the text which will later be used to specify the marked parts. Save your input.

#### Viewing and marking ECGs and measuring HR during a pause @

Clicking on Pause (2), will stop the screen display. The recording continues to run and is displayed on a channel (5). The scroll bar can be used to view the current recording (6).

The tools Mark, Measure HR and Measure can be found in the "Mouse Function" area. By dragging the red cursor in the ECG (using the Mark mouse function), you can mark sections. A dialogue appears for specifying, printing and saving the marked part. With Continue **1** you return to the normal view.

#### Online ECG printing (printing ECG)

By clicking on the Print <sup>(B)</sup> button a screen page of the ECG is printed from the point of clicking. The printout contains 4.5 to 9 seconds of the ECG, depending on the display speed <sup>3</sup>.

3) Note on online printing: Under Examination, Resting ECG, Settings, Print, General, you can define in the "Online ECG print settings" area whether the ECG should be printed as it appears on the monitor or whether online printing should be carried out according to previously defined print settings for the analysed ECG

#### 3.4.2 **Recording of stress ECG**

The profile selection is opened. Select a stress profile **1**. The list contains predefined profiles for ergometers (with watt indications) and the treadmill. Then set the training device for the recording **2**.

The predefined stress profiles can be modified and adapted as required. By clicking on Save 3 you can save the modified stress profile under a new name. With the New 4 button you can create new profiles (types: Ergometer, Free, Treadmill). The values in the Alert area **S** can be set freely and have to be activated if required.

If you want to use the alert function for stress ECG, the alert values in the profile selection must be set correctly before clicking on Start. The alert values cannot be activated or changed later.

After selecting and configuring the stress profile, click on the Start 6 button to access the recording screen.

#### Additional information: Steady state option for ergometer profiles

By using Steady state 2, the stress profile can be manually controlled during recording. If Steady state 7 is selected, entries can no longer be made for the stage duration, increase and end in the input mask. During the recording, the profile will continue running unchanged until a manual change is made. Click on the Measurement button to define the end of a load level during the recording. The last ten seconds will be measured. Then set the load for the new load level.

			User				custo	med GmbH				? _ >	<			? _
			Patie	nt			Muste	rmann Franz			10.10	.1960 (57 Y.)			10.10.	.1960 (
			Exam	nination			Stress	ECG					·			
Deskerel				2534-66	1			-		2	6003			2	. 6003	-
PTOLOCOI				25Wall +2	25Wall 2	min		•		ecsooo, Br Eigometer,	, SPO2	····· ·	1 1	Essembler	, 3102	
Resting phase	e			Stress phas	se					Recovery phase				Elgometer		
O End by us	ser			Stage		🔿 Rar	mp	7 O Stead	ly state	Number of stages	10					
Ouration		1:00	min	Stage dura	tion	2	min			Duration of 1st stage	2	min		Ergometer manuali	y	
Blood pre	ssure m	easureme	ent	Initial load		25	Watt			Load of 1st stage	25	Watt				
				Increase		25	Watt			Duration of next stage	2	min		Ireadmin manually		
Warming pha	ise			End ma	anually					Load of next stages	0	Watt		Blood pressure		
Duration		0:00	min	🔿 End at						Blood pressure 1st	stage			BP manually		
Load	0	0	Watt	Blood p	pressure m	1 2	min			Blood pressure me	2	min		BP Ergometer		
Blood pre	ssure													SPO2 measurement		
									Watt					SPO2 manually	_	_
									mate	Alert				Alert	460	
									135	Maximum HR	163	BPM		Maximum HR	163	BPM
									120	5 Maximum (systo	200	mmHg		Maximum (systo	. 200	mmF
									105	ST deviation	0.20	mV		ST deviation	0.20	mV
									90	Arrhythmia				Arrhythmia		
									75	Number VPB ma	3	BPM		Number VPB ma	. 3	BPM
									60	Patient		_		Patient		
									45	Age	57	Years		Age	57	Yea
									30	Height	175	cm		Height	175	cm
									15	Weight	70.0	kg		Weight	70.0	kg
										Target load	153	Watt		Target load	153	Wat
2		4	6	8	10	12		min		Target heart rate	143	BPM		Target heart rate	143	BPM
Start		Save				Delete			•		End				End	
								····· •								



10.10.1960 (57 Y

163 BPM 200 mmHg 0.20 m٧ BPM

If several ECG devices are connected to the workplace, the "Select ECG Device" (3) dialogue will be displayed. Select the ECG device (9) and click on Confirm (10).

If only one ECG device is connected, this step can be omitted.

Select ECG Device	8	7-593	
cardio 300 USB (s cardio 300 BT (Sf	5N: (9)300) N: UBE00300)	^	
Add	Confirm 10 Canc	el	

#### Monitoring and electrode control

The patient's ECG signal will be shown on the display but not recorded (monitoring). Work steps before the start:

- Change the type of lead if necessary 1.
- Check if all electrodes are attached optimally.
   If there are red lines on the screen, the contact between the skin and electrode(s) is insufficient. The corresponding electrodes will need to be reattached. For ECG devices with USB connection, selecting Electrodes (2) provides you with a graphic representation of the signal quality.
- > Set the required filter (Options menu (B).
- Recommended settings for stress ECG with treadmill: Under Options activate: Mains Filter, Muscle Filter and Ergo-Filter<sup>1)</sup>.
   ECG display Precordial, 5 mm/mV and 25 mm/s (line above ECG <sup>(1)</sup>).
- Click on Start <sup>(b)</sup> or Enter to start the recording.



1) The Ergo-Filter is only needed when strong movement artefacts are to be expected, e.g. when using a treadmill.

Keyboard shortcuts:



**Options** menu

#### Display and control elements (view after Start)

- Setting options for ECG display
- **b** Buttons for controlling and editing the ECG recording
- S Display of heart rate and blood pressure, countdown of the current level
- Change the current load and increase for ergometer profiles<sup>2)</sup>
   or speed and slope for treadmill profiles
- Load profile (orange) with heart rate curve (blue)
- **6** Blood pressure curve (green)
- Set the ST point, display the summary complex (selection of channel via buttons at left in front of the ECG signal)
- Show ST trend curve, ST values and overview of results (Online arrhythmia detection); button in the area
   flash red when the alarm limits are exceeded

#### Note on manual blood pressure measurement

You are regularly requested to measure blood pressure. Enter the values in custo diagnostic. In the Options menu, click on Blood Pressure or press the F9 key and enter the values. Click on Confirm to apply your input. The method of entering lactate (F10), SPO2 (F11) and Borg values (F12) is the same.



2) Note on the settings for the load change: You can define by how much watts Load and Increase should change each time the arrow buttons are presed

pressed. This setting can be found under: Examination, Stress ECG, Settings, Menu/functions, ECG view in the "Manual load change" area.



#### **Resting phase**

The resting phase begins after Start. This phase proceeds according to the settings in the profile selection, it has a minimum duration of ten seconds.

#### Stress phase

The stress phase then begins<sup>3)</sup>. This phase proceeds according to the profile. Manual load changes can be made at any time. The Next Stage <sup>(6)</sup> button can be used to end the current load level and start the next load level.

#### Entering an unconfirmed report during recording

Open the context menu and select **Report**. Enter the unconfirmed report in the large text field **1**. To save your input, click on **Confirm 1**. By pressing **Cancel 1**, the unconfirmed report is closed without any changes being applied.

If you save your information with Confirm <sup>(D)</sup>, the unconfirmed report becomes a (preliminary) report, depending on the reporting rights of the current user. The evaluation is thus (pre-)confirmed. If the evaluation is not to be classified as (pre-) confirmed at this point, you can reset the report status when selecting End.



		Reporting		an an	14	2/8	6 1:	22
V1		New unconfirmed report		<b>*</b>	Blo	A Pressu	Increase 25 W	•ss 1 •
V2 -		Ø						Load HR
V3 ~						1	6 8	min
V4		Descript ECG	ERGO-Val	▼ Assessm.				BF
VE	Λ ~ - Λ	F5 LoadTime F6 Nom.Load	F9 Abort F10 AngPect		- HF	t deper	ndent	▲ \/5
v3 -		F7 HF F8 BP-MAX	F11 Cyanosis F12 Normotor	1		A		00
V6		Options	Confirm	18 Cancel 19		l.	$\wedge$	

#### 3) Note on treadmill profiles:

The treadmill can be stopped using the Stop button, e.g. if a lactate measurement should be conducted. The treadmill will be restarted by clicking on the button again.

<u>Always warn the patient before you</u> <u>stop or start the treadmill!</u> If the Unconfirmed report option is active in the Settings, custo diagnostic generates an automatic unconfirmed report which is displayed in the evaluation. This option is enabled by default and can be disabled under Stress ECG, Settings, Diagnostic, Reports.

All unconfirmed reports produced by the system should be considered as suggestions only. For diagnosis and therapy purposes it is essential that the results are checked and assessed by a qualified physician.

#### Recovery phase, ending the recording

The recovery phase can be started using the Recovery phase <sup>(1)</sup> button, e.g. when the Manual end option was selected in the profile selection or as a result of a premature termination. Define the end of the stress phase (immediately or at the end of the load stage <sup>(2)</sup>). The dialogue for entering the reason for termination <sup>(2)</sup> is then opened <sup>4)</sup>. The reason for termination can be displayed in the evaluation.

If the end of the stress phase is defined in the profile, the recovery phase starts automatically after the last load level has expired. The recovery phase proceeds according to the profile.

If you would like to end the ECG recording but the ECG signal should still be displayed on the screen, click on Stop 2. Otherwise, the recording will be automatically saved, measured and displayed as an evaluation by clicking on the End button (bottom right).



4) Text modules for entering the reason for termination

Use: The text modules are called using the keyboard (F5 to F12) in the "Reason for termination" dialogue or by clicking on the corresponding button.

Configuring text modules: Select Examination, Stress ECG, Settings, Diagnostic, Reason for End to configure text modules for entering a reason for termination. Enter a name for the text module in the "Description" field. This name will appear on the button for calling the text module in the "Reason for termination" dialogue. In the "Text module" field, enter the text which will later be displayed as the reason for termination. Save your input.



#### Editing options during the recording process

#### Automatic marking of the ECG

Clicking on the Marker 29 button will automatically mark the last six seconds of the recording. A dialogue appears for specifying 5), printing and saving the marked part.

#### Viewing and marking ECGs and measuring HR during a pause

Clicking on Pause 29, will stop the screen display. The recording continues to run and is displayed on a channel 20. The scroll bar can be used to view the current recording 22.

The tools Mark, Measure HR and Measure can be found in the "Mouse Function" area 28. By dragging the red cursor in the ECG (using the Mark function), you can mark sections. A dialogue appears for specifying, printing and saving the marked part 29. With Continue 39 you return to the normal view.

#### Online ECG printing (printing ECG)

By clicking on the Print 3 button a screen page of the ECG is printed from the point of clicking. The printout contains 4.5 to 9 seconds of the ECG, depending on the display speed <sup>6)</sup>.

#### Text modules for specifying 5) marked parts

Use: If a part is marked during the recording, the "Markings" dialogue appears. The marked parts will be automatically specified by pressing the corresponding key (e.g. F5) or a previously configured text module button.

#### Configuring text modules:

Select Examination, Stress ECG, Settings, Menu/Functions, Markings to configure text modules for specifying marked parts. Enter a name for the text module in the "Description" field. This name later appears on the button for calling the text module. In the "Text module" field, enter the text which will later be used to specify the marked parts. Save your input.

#### Note on online printing: Select Examination

Stress ECG, Settings, Print, General to define in the "Online ECG print settings" area whether the ECG should be printed as it appears on the monitor or if online printing should be carried out according to previously defined print settings for the analysed ECG.





#### 3.4.3 Extended ECG settings

➤ Changing the ECG colour scheme: The ECG colours are preset in custo diagnostic and can be changed under Examination, Settings, System, ECG colour. Click on Save to apply your changes. The ECG grid in custo diagnostic corresponds to normal ECG paper. The small boxes measure 1 \* 1 mm, the large boxes 5 \* 5 mm. To ensure the graph paper is correctly displayed on the screen, the screen diagonal of the monitor must be specified in the custo service centre. Contact your authorised custo med dealer.

Resting ECG, automatic ECG procedures: It is possible under Examination, Resting ECG, Settings, Menu/Functions, Workflow to define the procedures for automatic ECG recordings in the "Automatic ECG" area. For example, recording duration and procedures after recording. Click on Save to apply your changes.

Procedures for manual resting ECG recordings and stress ECG: It is possible under Examination, Resting ECG or Stress ECG, Settings, Menu/Functions, Workflow to set the procedures after recording and the display options in the evaluation in the "Workflow" area. Click on Save to apply your changes.

> Print settings for resting ECG: It is possible under Examination, Resting ECG, Settings, Print, Print pages to define the contents for different printouts. In the "Type of printout" list, select the desired entry, for example automatic printing (automatic ECG). In the "Printout" area, select the contents for the printout after an automatic ECG.

Important: This setting is only required if the automatic ECG print pages contains content other than the standard printout (see "Type of printout" Standard). Click on Save to apply your changes.

Print settings for stress ECG: It is possible under Examination, Stress ECG, Settings, Print, Print pages to define the contents for different printouts. In the "Type of printout" list, select the desired entry and compile the contents of the print-out. Important: This setting is only required if the selected printout (see "Type of printout") contains content other than the standard printout (see "Type of printout" Standard). Click on Save to apply your changes.

Maximum load: The maximum load is displayed in the evaluation and serves for comparison with the nominal load. The criteria for determination of the maximum load are set on the Examination, Stress ECG, Settings, Diagnostic, Calculation screen. For example, when undercutting a certain duration the corresponding examination stage can be excluded.

#### Resting ECG, rhythm strips

During resting ECG recordings, additional rhythm strips can be recorded. These are ECG sections of any duration during which the recording can be manually controlled. The "Rhythm strip" function can be activated under Examination, Resting ECG, Settings, Menu/Functions, Workflow **1**. Define whether the recording is to be stopped manually **2** or automatically after a specified duration **3**. Save **4** your input.

The recording of a rhythm strip is activated on the ECG interface by clicking the Rhythm <sup>(5)</sup> button. If the recording duration is not set, the recording of the rhythm strip can be stopped using the Stop <sup>(6)</sup> button.

In the evaluation, the available rhythm strips can be called and displayed via the menu on the top left  $\bigcirc$ .

It is possible to print the rhythm strips. Select the Examination, Resting ECG, Settings, Print, Print pages screen and activate Rhythm <sup>(3)</sup>. Save your input.



### 3.5 Working with the evaluation

#### 3.5.1 Opening an evaluation via the evaluation search

- > To open the evaluation search<sup>1)</sup> right-click on the Patient button ①.
- With factory settings, the Search screen 2 is displayed.
   Here you can search for evaluations with previously defined and saved search criteria, so called filter sets. Filter sets can be created on the Advanced search screen 3.
- Depending on the default setting of the system, a filter set is already active and the search results are displayed here full-screen as a list 4.
- > If no filter set is active, select a set **5**.
- Open an evaluation by double-clicking on the corresponding line or via the Show evaluation button 6.

#### Configuring the list of search results

- Right-click on the screen to open the context menu.
   There click on Select columns and set the required columns.
   Click on Confirm to apply your changes.
- Clicking on a column header sorts by this column and sorting within the column can be reversed.
- > The list can be printed and exported in various formats 2.

#### Renaming filter sets, deleting filter sets

- Right-click on the screen to open the context menu.
   There click on Rename filter set oder Delete filter sets.
- > Follow the instructions.

#### **Editing filter sets**

> Open the Advanced search screen 3, see next page.

	User	custo med GmbH	? _ ×
	Patient 1		
	Examination		•
Search	2 Advance	d search 3	
current filter set empty	5		
	_		
Examination Date	Patient name	Patient ID Pr	operties
4			
-			
Evaluations: 0			
Print list 🚺 Export	Update list	Reassign evaluati	on Show evaluation 6
			Cancel

1) The evaluation search can be configured in the custo diagnostic settings, see Examination, Settings, Database, Eval. search.

#### Advanced search, creating filter sets

- The Advanced search screen ③ is used for creating filter sets and for quick selection of search criteria (e.g. examination, properties, period) ④.
   The search is limited by setting certain search criteria.
- > The search results are displayed in the lower half of the screen as a list @.
- An evaluation can be opened by double clicking on the corresponding entry in the list of search results or with the Show evaluation ① button.
- The previously selected search criteria can be saved for later use as a filter set with the appropriate name. Enter the name into the input field and click on Save current search as set B.

#### Editing filter sets

- > The filter set to be edited must be selected, see line "Current filter set".
- > Adapt the search parameters (e.g. examination, properties, period).
- Click on the Save current search as set <sup>B</sup> button to overwrite the previous set.
- > If a new name is assigned beforehand, a new set will be created.

#### Configuring the list of search results

- Right-click on the screen to open the context menu.
   There click on Select columns and set the required columns.
   Click on Confirm to apply your changes.
- Clicking on a column header sorts by this column and sorting within the column can be reversed.
- Use the arrow button ③ at the bottom right of the list to increase or decrease the size of the list.
- The list can be printed and exported in various formats 6.





Reference between the end dialogue and the evaluation search

In order to make proper use of the evaluation search, the status of the evaluation must be set correctly in the end dialogue when you exit an evaluation.

Example: An evaluation can only be found in the evaluation search with the property confirmed "No" if the status "Evaluation confirmed" is NOT selected in the end dialogue.

#### 3.5.2 Opening an evaluation via the evaluation main menu

- > Open the examination main menu via Resting ECG or Stress ECG.
- > Click there on Show Evaluation 1.
- The patient search screen appears. In this screen select the patient whose evaluation you want to open. Enter the patient's name into the input fields in the search screen 2.
- Select the patient from the list below the input fields 3
   and confirm your selection by clicking on the <u>Select Patient</u> 4 button or by double-clicking on the name.
- A list containing all the patient's evaluations is then displayed.
   Select the desired evaluation from the list and open it
   by means of a double-click or via the Show Evaluation button.

	User			User	
	Patient			Patient	
	Examination			Examination	
New Resting ECG		Last name		Mustermann 2	
		First name			
Show Evaluation		Patient ID			
Show Comparison		Patient Group		All patients	
		Assignment	Physician	All physicians	
Assign Evaluation			Physician ID		
Settings		Last name	First name	Date of birth	Pat. ID
		Mustermann	Absoluta	10.10.1960	000000001
		Mustermann	Franz	3 10.10.1960	0000000002
	C				
					6 patien
		Select Patient	4	Edit Patient	
		New Patient			
Cancel		Cancel			

#### 3.5.3 Resting ECG evaluation

The ECG evaluation is divided into two main areas, ECG and Measurement. The ECG screen is preset as the start screen, the Measurement screen can be alternatively set as the start screen if required <sup>1)</sup>.

From the sub-screens of the two areas the main screen of the other area can be accessed at any time.

1) Setting the start screen: The start screen of the evaluation can be set under Examination, Resting ECG, Settings, Menu/Functions, Workflow in the "Menu/Functions, Show Evaluation" area.



#### Navigation in the resting ECG evaluation

The buttons for opening further evaluation screens are located at the bottom of the screen. The labelling of the buttons changes as soon as you switch to a different evaluation screen. The button that has been clicked always contains the name of the screen you just left.

Example: You click in the evaluation (view: ECG as the start screen) on the Measurement button. You go to the Measurement evaluation screen and the previously clicked Measurement button changes to ECG. By clicking on ECG, you return to the ECG view again.



#### Resting ECG evaluation, start screen ECG

- Setting options for ECG display
- Mouse functions for detailed viewing and measuring of the ECG signal (Zoom, Analysis, Measure HR, Measure, Marking)
- G Further evaluation screens
- **d** Button for printing the evaluation<sup>2)</sup>
- Button for closing the evaluation

If the Measurement view is set as the start screen, the same operating and navigation elements can be found there (a) to (a).

The evaluation start screen can be set under: Examination, Resting ECG, Settings, Menu/Functions, Workflow in the "Procedure control, Show Evaluation" area. 2) Settings for the Print button:

Under Examination, Resting ECG, Settings, Menu/Functions, Workflow, in the "Workflow, Print" area, you can display the expanded print menu by clicking on the Print button (standard setting) and define whether printing should be conducted automatically without further settings according to the standard print settings (= defined printout).

The standard print settings for resting ECG can be found under Examination, Resting ECG, Settings, Print, Print pages.

Click on Save to apply your input.



#### Resting ECG evaluation, Options menu<sup>3)</sup>

- Print menu for temporary changes to the print settings
- **b** Button for exporting the evaluation (e.g. as Excel, PDF, DICOM...)
- C Assign the evaluation to another patient, if necessary
- **d** Manual blood pressure input (F9 key)
- HR trend, display of events in the ECG (e.g. VES)
- **f** Tables and graphics for heart rate variability<sup>4)</sup>
- New analysis of ECG signal for resetting manual changes in the ECG, additions to the report remain available
- Automatic creation of a new report after manual changes have been made in the ECG recording
- Switch filters in the ECG on or off (options: Display ECG as saved, unfiltered or filtered ECG – mains filter, muscle filter)
- Show or hide additional content in the right half of the screen:
   e.g. enlarged pacemaker spikes, HR curve and ST trend,
   summary complexes and report or table of measured values
- Show or hide pacemaker spikes
- Graphic flattening of ECG signal
- The bar below the ECG shows the heart rate instead of RR intervals in milliseconds

 Note on Options menu The scope and content of the Options menu changes based on which evaluation screen you are currently on.

In the Options menu under Measurement, you can activate the display of ST values and set which marked parts should be displayed in the summary complexes.

4) Note on RR variability For the RR variability to be displayed, at least five minutes of ECG needs to be recorded!

5) Only with recordings including high resolution ECG (custo cardio 300 or 400).



#### 3.5.4 Stress ECG evaluation

The ECG evaluation is divided into two main areas, ECG and Measurement. The ECG screen is preset as the start screen, the Measurement screen can be alternatively set as the start screen if required <sup>1</sup>.

From the sub-screens of the two areas the main screen of the other area can be accessed at any time.

1) Setting the start screen: The evaluation start screen can be set under Examination, Stress ECG, Settings, Menu/Functions, Workflow in the "Workflow, Show Evaluation".



#### Navigation in the stress ECG evaluation

The buttons for opening further evaluation screens are located at the bottom of the screen. The labelling of the buttons changes as soon as you switch to a different evaluation screen. The button that has been clicked always contains the name of the screen you just left.

Example: You click in the evaluation (view: ECG as the start screen) on the Measurement button. You go to the Measurement evaluation screen and the previously clicked Measurement button changes to ECG. By clicking on ECG, you return to the ECG view again.



"ECG" evaluation screen

"Measurement" evaluation screen

#### Stress ECG, start screen ECG

- O Setting options for ECG display
- Mouse functions for detailed viewing and measuring of the ECG signal (Zoom, Analysis, Measure HR, Measure, Marking)
- C Load profile with heart rate and blood pressure curve
- Tabular display of PWC<sup>2</sup> (physical working capacity) and MET (metabolic equivalent), further information on PWC and MET available in the appendix
- e Further evaluation screens
- **f** Button for printing the evaluation<sup>3)</sup>
- 9 Button for closing the evaluation

If the Measurement view is set as the start screen, the same operating and navigation elements can be found there (a) to (9), without (c) and (d).

The evaluation start screen can be set under: Examination, Stress ECG, Settings, Menu/Functions, Workflow in the "Workflow, Show Evaluation" area. 2) Setting the PWC predicted values: The PWC predicted values are preset in custo diagnostic and can be changed under Examination, Stress ECG, Settings, Diagnostic, Reference values. Click on Save to apply your input.

3) Settings for the Print button: Under Examination, Stress ECG, Settings, Menu/Functions, Workflow, in the "Workflow, Print" area, you can define whether the extended print menu should be displayed when the Print button is clicked (standard setting) or whether printing should be conducted automatically according to the standard print settings (= defined printout) without any additional settings being made.

The standard print settings for stress ECG can be found under Examination, Stress ECG, Settings, Print, Printed pages.

Click on Save to apply your input.



#### Stress ECG, Options menu<sup>4)</sup>

- a Print menu for temporary changes to the print settings
- **b** Button for exporting the evaluation (e.g. as Excel, PDF, DICOM...)
- C Assign the evaluation to another patient, if necessary
- **O** Display of blood pressure (F9), lactate (F10), SPO2 (F11) + Borg values (F12)
- G Trend graphics, e.g. on load, HR, ST, speed, blood pressure, lactate, SPO2...
- **f** HR trend, display of events in the ECG (e.g. VES)
- **9** Tables and graphics for heart rate variability<sup>5)</sup>
- New analysis of ECG signal for resetting manual changes in the ECG, additions to the report remain available
- Switch filters in the ECG on or off (Options: Display ECG as saved, unfiltered or filtered ECG – mains filter, muscle filter, Ergo-filter)
- Show or hide content in the right half of the screen:
   Summary complexes and report or table of measured values (preset:
   Trend = load profile with HR and blood pressure curve, PWC and MET)
- **(**Show or hide pacemaker spikes
- Graphic flattening of ECG signal
- The bar below the ECG shows the heart rate instead of RR intervals in milliseconds



In the Options menu under Measurement, you can activate the display of ST values and set which marked parts should be displayed in the summary complexes.

5) Note on RR variability For the RR variability to be displayed, at least five minutes of ECG needs to be recorded!



#### 3.5.5 Confirming the evaluation

#### Unconfirmed report and report

To open the unconfirmed report, right-click on the evaluation interface. In the context menu, select Report. Enter your data in the text field **1**. If the Unconfirmed report or Interpretation option is selected in the system settings, an automatic system unconfirmed report is already present in the text field. If necessary, older reports can be displayed via the report history (collapsible list above the text input field). When you click on Confirm **2** your input is saved and the unconfirmed report becomes a (preliminary) report, depending on the report rights of the current user. If your (unconfirmed) report is not yet complete but you want to save it nevertheless without reaching the "Evaluation (pre)confirmed" status, the report status is reset upon ending (End)the evaluation.

#### Text modules - an aid for writing reports

Select Examination, Resting ECG or Stress ECG, Settings, Diagnostic, Reports to configure text modules for confirming an evaluation **3**. A total of four groups **4** can be stored with up to eight text modules **5**. The text modules are called in the unconfirmed report dialogue using the keyboard (F5 to F12) **6**.

A text module can be created from normal text as well as variables. When you use a text module in the unconfirmed report, the actual value from the evaluation is inserted in the report text instead of a variable. The structure of a variable is  $\{VARIABLE\}$  (e.g. heart rate, resting:  $\{HR_R\}$ ). The Shortcuts for export values  $\heartsuit$  button provides you with a list containing all the variables. If the text modules should be shown in the unconfirmed report, make sure that the Enabled 3 option is activated.

Alternatively, the text modules can be shown in the unconfirmed report via Options, Texts on ②. You also have the option of entering a text, which will be automatically shown in each unconfirmed report ③. The text can be changed later in the unconfirmed report dialogue. Save your input.

				-			1	1	
00.01	•		0.05Hz - 125Hz	Resting ECG	•	Print	Menu/Functions	Export	Device
00.01	Unconfirmed Report		1008			Reports	Report Cont.	Calculation	Alarms
$\sim$	Current automatic unconfirmed repor	t by custo med GmbH, 29.01.2018 09							
~	normal heart rate sinus rhythm			Text modules for	report	3			Reporting
	normal type								Automatic rep
~	possible left ventricular hypertrophy			Category 4	- 1		Name Fre	eqRhyt	or interpretati
~	0		m	Function key 5	▼ F	5 🔺	Name no	rmRep	report status
			-1-	Text module	SR, R norm	esting HF: {HF al AV-conductio	_R}. Normal P-configura n (PQ {PQ_ZEIT}). Axis: 2/RAD/AT/LT/LAD_OPS	tion,	configure repo
~			-		norm (QTc:	al range ({QRS {QTC_ZEIT}).	ZEIT}). Regular repolar No irregularities.	istation	Measurement
~~~	FreqRhyt Axis	ECG Assessment	m		Prede	fined text mod	ules Shortcuts for e	xport values	
~	F5 normRep	F9 Rest-HF		Export elements	🔳 inc	lude units	V		Measurement
~	F6 normSR	F10	6						
	F7 Tachy	F11		Content of report	dialogu	9			User-defined
	F8 Brady	F12							
~~~	Options 9	Confirm 2 Cancel		Text modules		Hidden			0
90		90 90 90 90		Question	0	-nabled Hidden			
	, ,		~		•	Enabled			
	Measurement ECG Overview	Options	Print	Approval proc	ess				
		options							
				Save					

#### 3.5.6 Ending the evaluation

Click on End (bottom right) in the evaluation. The End dialogue opens. This is where the Status of Evaluation ① is defined <sup>1)</sup>.

- Evaluation preconfirmed: active if a user with the reporting right "Preconfirm evaluations" has confirmed the unconfirmed report.
- Confirmed: active if a user with the reporting right "Confirm evaluations" has confirmed the unconfirmed report. The "confirmed" status can be reset if required.
- **4 Printed**: indicates if the evaluation has been printed.
- Indelible: can be selected after reporting has been completed.
   The evaluation can now only be viewed and can no longer be changed.

Click on Confirm 6 to close the evaluation.

#### 3.5.7 Optional: Reporting with approval process

If custo diagnostic is used with approval process, then authorised persons with the corresponding user rights can save pre-reports of other persons as a report, without having to close the evaluation which was opened previously (shortened workflow) or enter pre-reports/reports directly if the evaluation was created by a person without reporting rights.

The approval process is visible in the unconfirmed report dialogue ① of an evaluation. The user or reporter can be changed there (User name ②, Password ③, Enter). During the logon process, the user rights of the respective user are checked and the software interface is adapted accordingly ③. The reporting is documented in the evaluation information ⑤ (context menu).

The approval process must be activated user and project-related in the settings and in the custo service centre. The user rights must be set according to the workflow <sup>2</sup>). Contact your authorised custo med dealer or custo med.

 Note: Pre-reporting physicians must have the user right Preconfirm evaluations, reporting physicians must have the user rights Confirm evaluations and Change reports of other users.

Muster

Sent via d

Received

imported

175cm / 70.0k

Patient

Examina

6

Assigned physician of patient: Patient height / weight / sex / age:

6

custo med GmbH

custo med GmbH

printed
indelible

Evaluation confirmed

 Activity
 Date
 User
 Worksta

 Status change 29.01.2018 09:15:49 custo med GmbH CDXM3N
 Woolfied
 29.01.2018 09:02:39 custo med GmbH CDXM3N

 Created
 29.01.2018 09:02:39 custo med GmbH CDXM3N
 Control of the transmission of transmission of the transmission of transmissicon of transmission of transmission of transmissico

Evaluation pre-confirmed Export to

Evaluation inform

Created by:

Preconfirmed by: Confirmed by:

Evaluation flag:

-	12	Channel 🔺 m	m/mV • 10	▲ mm/s ▼	50 🔺
J2 00 03 00					
kankanka	U	Inconfirmed Report	1	_	
		Current automa	tic unconfirmed rep	ort by custo med G	mbH, 29.01.20
	1	normal heart rat	e sinus rhythm		
tus of Evaluation		possible left ver	tricular hypertroph	у	
as of Evaluation	-				
j evaluation precommed					
confirmed					
printed	_				
indelible					
	-	FreqRhyt	Axis	ECG	Assessm
6 Confirm	~	F5 normRep		F9 Rest-HF	
Commit		F6 normSR		F10	
	-	F7 Tachy		F11	
		F8 Brady		F12	
i Ki	1-	Reporter	custo med GmbH	2 🔺	6
from from the	4	User rights: W	rite evaluation report, Pre	-confirm evaluations, Char	nge reports of other
	~	Ontions	•	Confirm	Cancel
		options		Commit	Cancer
90 90 90 90	7	5 75	75	75	75
90   90   90					
90 90 90 90					

1) The assignment of properties (status of evaluation) in the End dialogue makes it easier to find evaluations in the evaluation search.

#### 3.6 Appendix

#### 3.6.1 Values and formulas in the ECG evaluation

#### > PWC (physical working capacity)

The PWC values indicates the physical ability of a patient at a specific heart rate. The PWC value is specified in watt/kg (body weight). In custo diagnostic, the PWC value is determined for a heart rate of 130, 150 and 170. To determine the PWC value, the patient must have a heart rate of +/- 10 to the reference rate (130, 150 and 170). If the heart rate has not been precisely achieved, the PWC value will be calculated using interpolation or extrapolation. Example: If a patient who weighs 100 kg reaches a heart rate of 170 to 200, the PWC value will be calculated as follows:

#### PWC170 = 200 W : 100 kg = 2 W/kg

The PWC predicted values are preset in custo diagnostic and can be changed under Examination, Stress ECG, Settings, Diagnostic, Reference values. Click on Save to apply your changes.

#### MET (metabolic equivalent)

The metabolic equivalent is used to determine the expenditure of energy during the maximum load. In custo diagnostic, the metabolic equivalent is calculated as follows:

Treadmill ergometry:

- v = max. speed in miles per hour
- m = slope in %

#### MET = 1 + (v \* 26.8 \* (0.1 + m \* 0.018)) : 3.5

Bicycle ergometry:

- L = max. load in watts
- W = weight in kg

MET = 1 + (12 \* L) : (3.5 \* G)

#### Calculation of QTc duration

Formula according to Bazett:

$$QTc-Duration = QT * \sqrt{\frac{HR}{60}}$$
$$QTc-Duration = QT * \sqrt[3]{\frac{HR}{60}}$$

Formula according to Fridericia:

#### Calculation of the target load

custo diagnostic offers two different calculation options for calculating the target load at maximum workload: "standard" and "according to Prof. Froelicher". The settings can be found under Examination, Stress ECG, Settings, Diagnostic, Reference values.

Standard formula Male, under 30 years of age Female, under 30 years of age: Male, over 30 years of age: Female, over 30 years of age:

Target load = 3 \* weight Target load = 2.5 \* weight Target load = 3 \* weight \* ((130 - age) : 100) Target load = 2.5 \* weight \* ((130 - age) : 100)

Source: Rost, R. & Hollmann, W. (1982): Belastungsuntersuchungen in der Praxis Georg Thieme Verlag, Stuttgart, New York. 164 p.

#### Formula according to Prof. Froelicher

Female:

3.933 + (86.641 \* body surface area) - (0.015 \* age) - (0.346 \* body surface area \* age)

#### Male:

6.773 + (136.141 \* body surface area) - (0.064 \* age) - (0.916 \* body surface area \* age)

Comment: The body surface area is calculated according to the formula of DuBois & DuBois: BSA = 0.007184 + height [cm]0.725 • weight [kg] 0.425 Source: DuBois, D. & DuBois, E.F. (1916): A formula to estimate the approximate surface area if height and weight be known. Arch Intern Med, 17: 663

#### BORG values for stress ECG

When performing a stress ECG, it is possible to enter BORG values during recording. BORG values are used to evaluate the subjective perceived exertion and were established by the Swedish physiologist Gunnar Borg in the Borg scale named after him. Classification is carried out either by the physician or by the patients themselves.

#### Simplified scale of subjective perceived exertion:

06	very very light	13	somewhat hard
07	very very light	14	somewhat hard
08	very very light	15	hard
09	very light	16	hard
10	very light	17	very hard
11	light	18	very hard
12	light	19	very, very hard
		20	too hard, no longer possible

Simplified BORG scale: Source: http://www.uni-bielefeld.de (Christian Stallmann)

#### 3.6.2 Keyboard navigation and shortcuts in custo diagnostic

Use the quick links in the main navigation, the keyboard navigation and the keyboard shortcuts to enable fast and convenient working.

#### Quick links in the main navigation

User	custo med GmbH	0	? _	×
Patient 🚺		22		
Examination		88		-

#### LEFT-CLICK

2 Call last patient

**1** Change user's password

**B** Examination main menu

#### **RIGHT-CLICK**

- Evaluation search
- 2 Call last patient
- 3 Most recently opened evaluation

User	custo med GmbH 🛛  🕘	? _ ×
Patient	Mustermann Franz 🏮 4	10.10.1960 (57 Y.)
Examination	Resting ECG 65	•

#### LEFT-CLICK

#### **RIGHT-CLICK**

- **4** Change user's password
- **6** Patient master data
- Menu for the current examination

## 4 All evaluations for the patient5 Most recently opened evaluation

for this examination

#### **Keyboard navigation**

When you press the Alt key, the initial letter of all the buttons on a screen page is underlined. Pressing an initial letter again triggers the corresponding button.

	<u>U</u> ser	<u>c</u> usto med GmbH	? _ ×
	Patient		
	Examination		-
Holter			
Аврм			
Resting ECG			
<u>S</u> tress ECG			
Cardiopulmonary Exercise Testing			
Spirometry			
Telemetry			
<u>C</u> ardiac Rehab			
Mobile Cardiac Telemetry			
Blood glucose			
Ecg Streaming			
Task Manager			
Worklist		- custo di	

#### Generally valid keyboard shortcuts

- Esc End, cancel
- Enter Confirm, continue
- Ctrl I Program information
- Ctrl H User main menu
- Ctrl P Patient main menu
- Ctrl U Examination main menu
- Ctrl S Patient master data for the selected patient
- Ctrl A All examinations belonging to the selected patient
- Ctrl G List of the most recently opened evaluations (same as clicking on the arrow button at top right)
- Ctrl F List of the most recently opened evaluations
- Ctrl L Evaluation search
- Ctrl W Waiting room list
- Ctrl Q Device list
- Ctrl M Switch to Metasoft
- Ctrl F1 Create system report, service e-mail

#### Resting ECG keyboard shortcuts during the recording

- Enter Start the recording
  - Esc End the recording
    - + Increase amplitude
    - Decrease amplitude
  - F9 Blood pressure input dialogue

#### Stress ECG keyboard shortcuts during the recording

- Enter Start the recording
  - Esc End the recording
    - + Increase amplitude
    - Decrease amplitude
  - F7 Starts an additional blood pressure measurement
  - F8 Creates a new load level for steady state profiles
  - F9 Blood pressure input dialogue
  - F10 Lactate input dialogue
  - F11 SPO2 input dialogue
- F12 BORG input dialogue
- Arrow-keys right/left: increase/decrease the load rise (bicycle) or the slope (treadmill)
- Arrow-keys top/down: increase/decrease the load (bicycle) or the speed (treadmill)

#### Generally valid keyboard shortcuts in an open evaluation

- Ctrl N Unconfirmed report input dialogue
- Ctrl K Medication input dialogue
- Ctrl R Call comparison
- Ctrl T Call trend
- Ctrl D Call print dialogue
- Ctrl E Call settings
- Ctrl O Call options menu

## Manufacturer's contact details:

custo med GmbH Maria-Merian-Str. 6 85521 Ottobrunn Germany

 Phone:
 +49 (0) 89 710 98 - 00

 Fax:
 +49 (0) 89 710 98 - 10

 E-mail:
 info@customed.de

 Internet:
 www.customed.de

