



Operating Manual

Spirometry


1 Safety

2 Hardware

3 Software

4 Hygiene

Part 3: custo diagnostic software for custo spiro mobile



Operating characteristics:
custo diagnostic 5.0
and higher
for Windows®

MSW 0006 – DK 1688
Version 002 – 16/01/2019

CE 0123

 **custo·med**
EXCELLENCE IN DIAGNOSTICS



© 2019 custo med GmbH

This Operating Manual may not be copied in its entirety or in part, duplicated in any form or by any means or translated into another language without the prior written consent of custo med GmbH.

The manufacturer reserves the right to change the information in this Operating Manual without prior notice. The current version can be downloaded from our website: www.customed.de.

CAUTION:

This Operating Manual is part of a modular system, consisting of four parts. All four parts must be downloaded from the Internet or from a CD to ensure the Operating Manual is complete.



Operating Manual

Spirometry

1 Safety

2 Hardware

3 Software

4 Hygiene

Part 3: custo diagnostic software for custo spiro mobile

Table of contents

3.1	Symbols used in this Operating Manual	4
3.2	custo diagnostic program structure	5
3.3	Connecting custo med devices to the PC	6
3.3.1	Connecting custo spiro mobile to the PC	6
3.3.2	Calibrating custo spiro mobile	7
3.4	Performing the spirometry measurement with custo diagnostic	
3.4.1	Reference measurement	8
3.4.2	Follow-up measurements: Spasmolysis and provocation	15
3.4.3	Unconfirmed report	17
3.4.4	Printing the measurement	18
3.5	Working with the evaluation	
3.5.1	Opening an evaluation	19
3.5.2	Structure of reference and spasmolysis measurements	21
3.5.3	Structure of provocation measurements	22
3.5.4	Navigation in the evaluation	23
3.5.5	Reference and spasmolysis evaluation	24
3.5.6	Provocation evaluation	25
3.5.7	Additional screens of a spirometry evaluation	26
3.5.8	Confirming the evaluation	29
3.5.9	Ending the evaluation	30
3.5.10	Optional: Reporting with approval process	30
3.6	Settings for the spirometry	
3.7	Error messages and remedies	36
3.8	Appendix	
3.8.1	Measured values in custo diagnostic	37
3.8.2	Calculation tables for predicted values	39
3.7.3	Keyboard navigation and short-cuts in custo diagnostic	52

3.1 Symbols used in this Operating Manual

ACTIONS THAT ARE PROHIBITED

or not allowed under any circumstances!



WARNING

used to indicate situations which, if not avoided, could result in personal injury or property damage



NOTE

provides important information which must be observed



TIP

contains practical information to assist you with your work



Words highlighted in colour indicate buttons or click paths to the corresponding program point, e.g. Examination, Spirometry

Words highlighted
in colour...

3.2 custo diagnostic program structure


The custo diagnostic program is divided into three areas – **User**, **Patient** and **Examination**. This structure ensures that you can always recognise who (which user) is carrying out what type of examination with whom (which patient). The main menus of each area can be reached by clicking on **User**, **Patient** or **Examination**.

In the **User** ① main menu, the user of the system can be selected. The administration of users takes place in the custo diagnostic service centre (creating users, allocating user rights, user-specific settings).

The **Patient** ② main menu is used for patient management. Its most important functions include **Search for patient**, **Add new patient** and **Search Evaluation**.

The **Examination** ③ main menu lists all of the examination types which are possible with custo diagnostic. Modules already purchased are active (black font), all other modules are inactive (light grey font).

This menu is also linked to the **Settings** ④ area. This area is for making cross-program, examination-related and user-specific settings.



3.3 Connecting and configuring the device

Prerequisite: custo diagnostic is installed on your PC and ready to use. The custo med devices and components may only be connected to the PC after custo diagnostic has been installed. The required device drivers are installed on the PC via the custo diagnostic standard setup or by specific selection during the custo diagnostic setup.





3.3.1 Connecting custo spiro mobile to the PC

- Connect custo spiro mobile to the PC. Driver installation takes place automatically. Wait until the installation is complete.
- Start custo diagnostic. Open the **Examination, Spirometry, Settings, Device, Device connection** screen ①.
- Activate the **spiro mobile** option ②.
- If working with custo spiro protect bacterial and viral filters, activate the **spiro protect** option. ③
- Enter the code for the custo spiro protect bacterial and viral filters ④. The five-digit code can be found on the packaging of the filters (refer to the Cat. No. line).
- Click on **Save** ⑤ to apply your input.
- Click on **End** ⑥ to close the settings screen.

3.3.2 Calibrating custo spiro mobile

- Before using the device for the first time, a calibration must be carried out. Observe the instructions provided in the hardware description.
- Connect custo spiro mobile to the PC.
Connect the devices together as illustrated ①.
- Start custo diagnostic and open the Examination, Spirometry, Calibration ② screen.
- Specify the volume of your pump ③.
- Enter your name in the performed by ④ field.
- Click on Start Calibration ⑤.
- Follow the instructions issued by the system ⑥.
- **After the start, pump rapidly from end stop to end stop. If your speed is not optimal, you will receive precise instructions.**
- As soon as the exact pump volume has been measured, the calibration process is automatically ended.
- The last calibrations are listed under "Documentation ..." ⑦.
- The Print calibration entries ⑧ button can be used to print a list of calibrations performed so far.
- The Delete calibration entries ⑨ button can be used to delete older entries from the list as required (delete entries that are older than ...).



3.4 Performing the spirometry measurement with custo diagnostic

Note on the procedure

The steps necessary to perform and evaluate a spirometry measurement in custo diagnostic are shown without a surgery IT system or HIS connection.



3.4.1 Reference measurement

Program start, calling the spirometry

- Start custo diagnostic and log in.
- Click on **Examination** ①, **Spirometry** ②, **New Spirometry** ③.

Selecting the patient

- Select a patient for the examination:
Enter the patient's name into the input fields in the search screen ④.
- Select the patient from the list.
Confirm your selection with **Select Patient** ⑤.
You can also select the patient by double-clicking on the name.

New patient

- If the patient does not yet exist in your database: Click on **New Patient** ⑥.
- Enter the patient data. The fields marked with an asterisk are mandatory.
In order to perform a spirometry measurement, you have to enter information on the patient's age, gender, height, weight and ethnicity¹⁾.
- **Save** the data; the patient is entered into the database.
- The settings for the reference measurement are then opened.

1) This information is used to determine the predicted value and to calculate the predicted values individually. If the required information has not been stored, you will be prompted to enter it later.

The first screenshot shows the main menu with 'Examination' (1) and 'Spirometry' (2) highlighted. The second screenshot shows the 'New Spirometry' (3) screen. The third screenshot shows the patient selection screen with a search bar (4) and a list of patients. The 'Select Patient' (5) and 'New Patient' (6) buttons are at the bottom.

Last name	First name	Date of birth	Pat. ID
Musterfrau	Martina	10.10.1978	0000000005
Mustermann	Absoluta	10.10.1960	0000000001
Mustermann	Franz	10.10.1960	0000000002

Reference measurement settings

- **Predicted value ①:** The default setting is GLI (Global Lung Initiative)²⁾.
The area of validity is displayed at the bottom of the screen ②.
If the patient data deviate from the area of validity, the system will propose a suitable predicted value. If no suitable predicted value is available, "none" is displayed in the "Predicted value" field. In this case, select a predicted value yourself from the menu. The determined predicted values are displayed in brackets during the measurement and may be incomplete.
- **Smoking habit ③:** Select the corresponding entry from the menu.
The "Smoking habit" option is included on the printout in the "Unconfirmed report" field.
- **Measurement type ④:** The first measurement taken from a patient is always a reference measurement. The measurement type can only be changed with follow-up measurements (spasmolysis or provocation).
- **Medication and Dosage (optional) ⑤:** You can specify here what medication was administered to the patient. For this function to be available, a selection of possible medications must have been previously stored in custo diagnostic. To do so, open the **Examination, Spirometry, Settings, Diagnostic, Drugs** screen. Enter the medication in the "Drugs" area ⑥ and click on **Add ⑦**. The procedure is identical in the "Dosage" area. **Save** your input.
- **Test with custo spiro protect ⑧:** To obtain precise measurements, indicate if you are using bacterial and viral filters.
- **Environment Data ⑨:** Set the environment data correctly (weather station, Internet). **Confirm** your input.

2) The preset predicted value can be changed in the custo diagnostic settings, see chapter 3.6 Settings for the spirometry.

GLI is characterised by the following features:

- The underlying population for the predicted values is much more extensive and broad than for other predicted values.
- The equation for determining the average predicted values is more accurate than for other predicted values.
- There is no noticeable transition from childhood to adulthood.

② GLI: Adults: 18 to 95 Years, 145 to 200 cm, ethnicity: Caucasian, Black, South East Asian, North East Asian

Overview of the measurement interface

- The predicted value that has been selected for the measurement series, in this case GLI
- Orientation aid (only for GLI, otherwise predicted value curve) constructed from FVC and FEF25-75, with the display of the predicted value range (green bars) formed from the predicted value and LLN³⁾
- Display of the results for FEV1, FVC and FEV1/FVC in a bar diagram; arrows mark the respective result after the measurement. Values within the green ranges can be considered as acceptable. Values located in the grey areas of the bars are considered pathological, divided into light, moderate and severe.
- During the measurement: Inclination sensor to check the posture of the patient⁴⁾
After the measurement: Miniature views of the performed measurements
- Table of measured values with predicted values, measured values obtained, Z-score⁵⁾ and measured value deviations in percentage from the predicted values
- Instructions for performing the breathing manoeuvre
- Display of the time-volume chart in real-time
- Settings for the measurement (change predicted value, animation on/off)
- Starting or stopping the measurement
- Ending the measurement, closing the measurement interface



Note on the functionality:

The Z-score, LLN, report assessment and the explanation according to clinical and occupational criteria are only available for measurements with the GLI predicted value.

3) LLN (Lower Limit of Normal) is the lower limit value, used to assess "normal" or "pathological". LLN corresponds to the 5% percentile of a healthy population. This means that if a measured value is below the 5% percentile, there is a 95% probability that a pathological finding exists or a patient with the corresponding value is healthy in only 5% of the cases. The green bars in the area of the measurement curve or in the coordinate system – orientation aid (b) – are formed from the predicted value (upper edge) and LLN (lower edge). Measurement curves above or within the green bars can be considered acceptable. The same applies to all other green bars in the spirometry surface. Measured values \geq LLN, i.e. within the green range, are considered acceptable.

4) The "Inclination sensor" function is part of the professional software and not included in the standard scope of supply. If no inclination sensor is used, this area remains empty during the first measurement. The miniature views of the measurements that have already been performed are then displayed.

5) The Z-score indicates by how many standard deviations a certain measured value deviates from the average predicted value. For example, $Z = 0$ corresponds exactly to the average predicted value and $Z = -2$ means that the measured value is two standard deviations below the average predicted value. A specific percentile can always be assigned to each Z-score. A Z-score of -1.645 corresponds to the 5% percentile (LLN). If the Z-score is greater than or equal to -1.645 , the measured value is not in the pathological range. The Z-score for the corresponding measured value is marked with a green square in the table of measured values. If the Z-score is smaller than -1.645 , the value is marked with an orange-coloured square (see the guideline on spirometry. Pulmonology. 2015; 69: 146-163).






Optional: Spirometry measurement with animation for children⁶⁾

To perform a measurement with animation for children, click on **Options**, **Animation** **k**. The procedure for performing a spirometry measurement is now demonstrated by means of animated drawings: the little dragon tries to spit fire in order to light a candle **l**. With his/her breathing, the patient supports the dragon to light the candle⁷⁾. The procedure with animation can be shown and explained before the start of the measurement via **Options**, **Instruction** **m**.

The animation should run on a separate screen for the patient. To set an additional screen, open the **Examination**, **Spirometry**, **Settings**, **Menu/Functions**, **Animation** **n** screen and select the **Own window** **o** option in the "Animation" area. On this screen it is also possible to set for which age groups the animation should be automatically displayed **p**. **Save** your input (left bottom button).

6) The "Animation for children" function is part of the "professional" software and not included in the standard scope of supply.

7) Under **Options**, **Setting**, **Animation** you can define the exhalation level (PEF) at which the candle is lit. The lower the percentage, the easier it is to light the candle. Click on **Confirm** to apply the changes.



Performing a reference measurement


- The patient must wear a nose clip.
- Click on **Start** ①.
- Next, the patient places the mouthpiece of the device into his/her mouth.
The mouthpiece must be held firmly between the patient's lips.
Provide clear instructions on how to perform the breathing manoeuvre⁸⁾.
The maximum recording time is two minutes.
- Remove the device from the patient after the breathing manoeuvre.
- To cancel an on-going measurement, e.g. if the device is handled incorrectly or the breathing manoeuvre is performed wrongly, click on the **Stop** button.
- Measurements that have been performed are displayed in the form of **miniature views** ② (several if the **Repeat** function was used).
- The most recent measurement is displayed in the coordinate system ③.
Other measurements can be displayed by clicking on the corresponding **miniature views** ②.

8) Before the breathing manoeuvre takes place, the patient must take a few resting breaths. The number of resting breaths can be changed in the settings.
To do so, open the **Examination, Spirometry, Settings, Diagnostic, Parameters** screen. The setting can be found in the "Breathing manoeuvre" area. **Save** your input.

Additional feature inclination sensor

The inclination sensor ① is displayed after pressing the **Start** button. This feature is used to check the posture of the patient during the measurement. An upright posture improves the quality of the measurement and allows a more precise identification of the measured values, in particular, FEV1.

The white arrows of the inclination sensor ② should be located in the middle, green area of the display. If the patient leans too far to the front or rear, the arrows will move into the red or yellow area, which may affect the determination of the measured values. The miniature views of the measurements that have been performed have a green, yellow or red background ③ in order to provide information about the patient's posture during the measurement. Green: upright posture (-10° to $+10^\circ$), yellow: inclined too far back, red: inclined too far forward.



Repeat measurement

- Click on **Repeat** 4 in order to perform another reference measurement. Up to six repetitions are possible.
- If more repeat measurements are required (more than six), measurements in the miniature views area will need to be deleted.
- To do so, click with the left mouse button on the measurement that you would like to delete 5, open the **context menu** by right-clicking and click on **Delete Meas.** 6.

Defining the best measurement

- The best measurement in a series of measurements is determined by the software⁹⁾ and marked in bold in the same colour as the measurement curve 7.
- To define another measurement as the best measurement, click with the left mouse button on the miniature view of the desired measurement, open the **context menu** by right-clicking and click on **Best Meas.** 8.

9) The best measurement is determined from the sum of FVC and FEV1. The best measurement can also be determined using other values. This setting can be found on the **Examination, Spirometry, Settings, Diagnostic, Parameter** screen in the "Identification best value" area.



Documenting patient cooperation

- Open the **context menu** by right-clicking, select **Cooperation** 9 and evaluate the patient's cooperation. Click on **Confirm** to transmit your data into the unconfirmed report.

Ending and closing the measurement

- With the default settings, the measurements are checked for reproducibility. If there are two reproducible measurements¹⁰⁾, a corresponding message 10 appears and the measurement can be finished.
- Starting a **follow-up measurement**: If spasmolysis or provocation is to be carried out immediately after the reference measurement, click on **Follow-up measurement** 11.
- Click on **End** (bottom right) and then on **Confirm** 12 in the End dialogue.

10) This function can be switched on and off in the settings. In addition, the required number of reproducible measurements can be changed and the criteria for reproducibility adapted. These settings can be found on the **Examination, Spirometry, Settings, Menu/Functions, Menu/Functions** screen in the "Flow Control" area at the very bottom. **Save** your input.









Options during the reference measurement

- During the reference measurement, the following editing and reporting functions are available in the **Options a** menu:
 - b** Print...: Print menu for compiling a printout
 - c** Change the Predicted value
 - d** Information on the spirometric lung age¹¹⁾,
 - e** Autom. Report¹²⁾,
 - f** Explain Report: Table with measured values, predicted values, limit values and the Z-score. Explanation of the assessment criteria underlying the assessment functions clinical, occupational and COPD-GOLD¹³⁾.
 - g** Report evaluation on/off: Assessment of the measurement results in a bar diagram below the measurement curve
 - h** Progress: Superimposition of the measurement curves to check plausibility

11) The spirometric lung age is determined using the FEV1 value, if FEV1 is calculated depending on age (not for all predicted values). The spirometric lung age is calculated based on the deviation from the predicted value. The spirometric lung age can be displayed via **Options, Lung age**. On the **Examination, Spirometry, Settings, Diagnostic, Parameter** screen you can define which predicted value should be used to calculate the spirometric lung age.

12) The following types of unconfirmed reports can be selected: Standard (70% rule for FEV1/FVC and 80% rule for IVC and FVC), COPD-GOLD, clinical or occupational evaluation according to GLI.

13) Z-Score, LLN and assessment criteria clinical and occupational only for GLI.



3.4.2 Follow-up measurements: Spasmolysis and provocation

Spasmolysis and provocation are referred to as sequential measurement. These measurement types can only be carried out after a reference measurement.

Calling the follow-up measurement

- Start the follow-up measurement via **Examination, Spirometry, Follow-up measurement** ①.
- A list will open including all patients for whom a reference measurement has already been carried out on the current day (= Spirometry group ¹⁾) ②.
- Select the patient from the list ③ and confirm your selection with **Select Patient** ④ or double-click on the name.

1) Note on the spirometry group:
A patient is deleted from the spirometry group if a follow-up measurement has been saved, if a new measurement has been generated (via **Examination, Spirometry, Measurement**) or on the next day.

Sequential measurement settings

- Define the settings for the follow-up measurement.
Check and, if necessary, change the data.
- Select the measurement type – spasmolysis or provocation ⑤.
- You can only enter data for "Medication" ⑥ if a medication list has been previously created in the settings for spirometry (**Examination, Spirometry, Settings, Diagnostic, Drugs**).
You can also carry out a follow-up measurement without indicating the medicine administered.
- Click on **Confirm** ⑦ to continue.

Note on follow-up measurements

The setting options and control elements for follow-up measurements are the same as for a reference measurement, e.g. **Start**, **Repeat** and **Best Meas..**

**Performing a spasmolysis**

- The patient must wear a nose clip.
- Click on **Start**.
- Next, the patient places the mouthpiece of the device into his/her mouth. The mouthpiece must be held firmly between the patient's lips.
Provide clear instructions on how to perform the breathing manoeuvre.
The maximum recording time is two minutes.
- Remove the device from the patient after the breathing manoeuvre.
- **Repeat** 8 the measurement if necessary. Measurements that have been performed are displayed in the form of **miniature views** 9 at the bottom left.
- The results of the spasmolysis (orange) are shown together with the results of the reference measurement (blue) for direct comparison 10.
- To close the measurement interface, click on **End** (bottom right).



Performing a provocation test²⁾

- The patient must wear a nose clip.
- Click on **Start**.
- Next, the patient places the mouthpiece of the device into his/her mouth. The mouthpiece must be held firmly between the patient's lips.
Provide clear instructions on how to perform the breathing manoeuvre.
The maximum recording time is two minutes.
- Remove the device from the patient after the breathing manoeuvre.
- **Repeat** the measurement if necessary.
- The results of the provocation (green) are shown together with the results of the reference measurement (blue) for direct comparison 11.
- Under **Options**, **Provocation test** the PD20 provocation dosage is displayed³⁾.
- The next measurement is triggered via the **Follow-on measurement** 12 button (or later via **Examination, Spirometry, Follow-on measurement**).
- To close the measurement interface, click on **End** 13.

2) Procedure for performing a provocation test: A maximum of eight provocation measurements are followed by one dilation measurement and up to eight control measurements.
The different measurement types for a provocation test are highlighted in colour in custo diagnostic to enable them to be distinguished more easily:

Reference measurement — blue
Provocation — green
Dilation — orange-brown
Control measurement orange-brown

3) PD20 provocation dosage: Medication dosage for the 20 percent drop of FEV1 in a provocation measurement compared to the initial value in the reference measurement.



3.4.3 Unconfirmed report

To open the **Unconfirmed Report** ❶, right-click on the measurement interface and select **Report** in the context menu.

If the **Unconfirmed report** option is enabled in the system settings, the unconfirmed report dialogue will already contain an automatic system report ❷. This option is enabled by default and can be disabled under **Spirometry, Settings, Diagnostic, Reports**. You can modify and supplement the text in the unconfirmed report dialogue. To save your input, click on **Confirm** ❸. By pressing **Cancel** ❹, the unconfirmed report is closed without any changes being applied.

If you save your input with **Confirm** ❸, the unconfirmed report becomes a (preliminary) report, depending on the reporting rights of the current user. The evaluation is thus (pre-)confirmed. If the evaluation is not to be classified as (pre-)confirmed at this point, you can reset the report status when selecting **End** in the End dialogue.

All unconfirmed reports produced by the system should be considered as suggestion only. For diagnosis and therapy purposes it is essential that the results are checked and assessed by a qualified physician.



Additional automatic reports

By selecting Options, Autom. Report, you can add the following evaluations to the report ❺:

- Standard, according to the 70% rule for FEV1/FVC and 80% rule for IVC and FVC
- COPD-GOLD¹⁾, information about the presence of chronic obstructive pulmonary disease and its severity
- Clinical evaluation according to GLI
- Occupational evaluation according to GLI

1) custo diagnostic can be set so that the COPD-GOLD report is automatically included in the unconfirmed report. To do so, open the **Examination, Spirometry, Settings, Diagnostic, Autom. Report** screen and activate the **COPD-GOLD** option. Specify whether the COPD-GOLD report should be included after a **reference measurement** or after a **spasmolysis** in the unconfirmed report. The other automatic reports can also be configured here. **Save** your input.

	GLI	Pred	Ref.
IVC	3.56	3.53	

Unconfirmed Report ❶

Current automatic unconfirmed report by custo med GmbH, 30.01.2018 1

Report clinical evaluation according GLI:
normal

Report Standard:
normal.

collabora

F5	Optimal	F9
F6	Understand	F10
F7	Pain	F11
F8	Cough	F12

Options

Confirm ❸ Cancel ❹

	GLI	Pred	Ref.
IVC	3.56	3.53	
FVC	3.56	3.38	
FEV1	2.93	2.88	
PEF (l/s)		6.36	

Selection automatic report ❺

Add new automatic report.

automatic report COPD - GOLD

By adding a new automatic report, the existing confirmed report will change. Please check the report afterwards.

Confirm Cancel

4
3
2
1

4 8

3.4.4 Printing the measurement

► Print with system settings:

In the measurement interface, click on the **Print**¹⁾ **1** button.

► Print with temporarily changed print settings:



If you do not want to print with the system settings, open the **Options, Print...** screen in the measurement interface **2** and carry out your print settings there. Changes in this print menu only affect the current print.

Start printing by clicking on the **Print**³ **3** button.

► Print with the Job Manager:

The print orders will be collected in the Job Manager from where they can be printed collectively at a later point in time. In order to store print jobs in the Job Manager, open the **Options, Print...** **2** screen and click on **Print Task** **4**. The existing print jobs are started on the **Examination, Job Manager** using the **Execute/Execute all** button.

1) The system settings for the printed pages of a spirometry measurement can be found in custo diagnostic on the **Examination, Spirometry, Settings, Print, Printed pages** screen. On the **Examination, Spirometry, Settings, Print, General** screen you can define in the "Print sequence control" area which printed pages are printed when the **Print** button is pressed (Current page, Pre-set pages or Total Printout PA).



3.5 Working with the evaluation

3.5.1 Opening an evaluation

custo diagnostic offers different options to open an evaluation, e.g. via the evaluation search or the main menu of the respective examination.

➤ Evaluation search

Right-click on the **Patient** ① button. This opens the evaluation search.

In the Examinations area, enter what type of evaluation you are searching for, e.g. **Spirometry** ②. In the Properties area ③ you can define more search criteria.

Example: If you set the "confirmed" property to **No**, you will receive a list of all the evaluations which have not yet been confirmed – a type of to-do-list.

To start the search, click on **Search evaluations** ④ or activate **Search automatically** ⑤. This option triggers an automatic search in your database whenever the search criteria are changed.

The right half of the screen displays a list of all the evaluations which correspond to the activated search criteria. To open the desired evaluation, select it from the list and click on the **Show Evaluation** ⑥ button or double-click on the evaluation.

If you want to use the same search criteria for the next search, activate the **Save selection** ⑦ option.

Exam.	Date	Patient	Prop.
Spirometr	30.01.2018	Musterfrau Martina	-----
Spirometr	30.01.2018	Musterfrau Martina	-----
Spirometr	20.11.2017	Musterfrau Martina	-----
Spirometr	16.11.2017	Musterfrau Martina	-----



Reference between the end dialogue and the evaluation search

In order to make proper use of the evaluation search, the status of the evaluation must be set correctly in the end dialogue when you exit an evaluation.

Example:

An evaluation can only be found in the evaluation search with the property confirmed "No" if the status "Evaluation confirmed" is NOT selected in the end dialogue.

► Examination main menu

Open the examination main menu via **Examination, Spirometry** and click there on **Show Evaluation** ①.

Last name	First name	Date of birth	Pat. ID
Musterfrau	Martina	10.10.1978	0000000005
Mustermann	Absoluta	10.10.1960	0000000001
Mustermann	Franz	10.10.1960	0000000002

The patient search mask appears. In this screen select the patient whose evaluation you want to open. Enter the patient's name into the input fields in the search screen ②.

Select the patient from the list below the input fields ③ and confirm your selection by clicking on the **Select Patient** ④ button. You can also select the patient by double-clicking on the corresponding name.


A list containing all the patient's evaluations is then displayed. Select the desired evaluation from the list ⑤ and open it by means of a double-click or via the **Show Evaluation** ⑥ button.

Examination	Date	Time	Status
Spirometry	30.01.2018	10:11	B-----
Spirometry	30.01.2018	10:09	-----
Spirometry	30.01.2018	09:53	-----
Spirometry	20.11.2017	11:22	B-----
Spirometry	20.11.2017	10:16	-----
Spirometry	16.11.2017	09:11	-----

3.5.2 Structure of reference and spasmolysis measurements

Overview of the reference and spasmolysis measurements:

Display of the best measurement(s) each with a measurement curve(s) and a table of measured values



Note on the functionality:

The Z-score, LLN, report assessment and the explanation according to clinical and occupational criteria are only available for measurements with the GLL predicted value.



1) All evaluations for a patient as a bar diagram with a table of measured values. In addition, the Z-score trend can be displayed for each measured value.

2) Indicates the probability of developing a disease and its severity.

3) Specified in years, calculated based on FEV1 depending on age, not possible for all predicted values.

4) The following types of unconfirmed reports can be selected and added to the report text: Standard (70% rule for FEV1/FVC and 80% rule for IVC and FVC), COPD-GOLD, clinical or occupational evaluation according to GLL.


5) Assessment criteria for clinical, occupational and COPD-GOLD reports

6) PD20 provocation dosage: Medication dosage for the 20 percent drop of FEV1 in a provocation measurement compared to the initial value.

3.5.3 Structure of provocation measurements

Overview of provocation

Bar diagram and table of measured values with reference, provocation, dilatation and control measurements



Note on the functionality:

The Z-score, LLN, report assessment and the explanation according to clinical and occupational criteria are only available for measurements with the GLL predicted value.



1) All evaluations for a patient as a bar diagram with a table of measured values. In addition, the Z-score trend can be displayed for each measured value.

2) Indicates the probability of developing a disease and its severity.

3) Specified in years, calculated based on FEV1 depending on age, not possible for all predicted values.


4) The following types of unconfirmed reports can be selected and added to the report text: Standard (70% rule for FEV1/FVC and 80% rule for IVC and FVC), COPD-GOLD, clinical or occupational evaluation according to GLL.

5) Assessment criteria for clinical, occupational and COPD-GOLD reports

6) PD20 provocation dosage: Medication dosage for the 20 percent drop of FEV1 in a provocation measurement compared to the initial value.

3.5.4 Navigation in the evaluation


The buttons for opening the various evaluation screens are located at the bottom of the screen. By pressing one of the buttons, e.g. **Comparison** ①, the comparison view is opened and the name of the button changes to **Evaluation** ② (the name of the previous screen). When you click on the **Evaluation** ② button you can return to the overview.



"Overview"
evaluation screen

"Comparison"
evaluation screen

3.5.5 Reference and spasmolysis evaluation



- a Flow-volume charts
(reference measurement: blue, spasmolysis: orange)
- b Predicted value, in this case GLI
- c Display of the results for FEV, FVC and FEV1/FVC in a bar diagram;
blue or orange arrows mark the respective result.
Values within the green ranges can be considered as acceptable¹⁾.
Values located in the grey areas of the bars are considered pathological,
divided into light, moderate and severe. The evaluation can be performed
according to the criteria clinical, occupational or COPD-GOLD,
refer to the options below the diagram
- d Table of measured values with predicted values, measured values, Z-score²⁾
and deviations in percentage; clicking on the column header
shows or hides the corresponding curve (a)
- e Volume-time curve
- f Environmental data input before the measurement
- g Comparison between the patient's current and previous evaluation
- h Comparison of individual measurements from a series of measurements
(plausibility check)
- i Reduced report assessment (c) and, if applicable, resistance results
- j Options menu: Print menu, Trend (all evaluations for the patient as a bar
diagram and Z-score trend of all evaluations),
Export evaluation (e.g. Excel, PDF, e-mail), Change predicted value,
Medication Overview, Miller Quadrant³⁾, Lung age⁴⁾, Autom. Report⁵⁾,
Explain Report⁶⁾, Reallocate evaluation
- k Print with system settings
- l End evaluation

1) The green ranges in the bar diagram are defined by the predicted value (upper limit/right end) and LLN - Lower Limit of Normal (lower limit/left end). LLN is the lower limit value used to assess "normal" or "pathological". LLN corresponds to the 5% percentile of a healthy population. This means that if a measured value is below the 5% percentile, there is a 95% probability that a pathological finding exists or a patient with the corresponding value is healthy in only 5% of the cases.

2) The Z-score indicates by how many standard deviations a certain measured value deviates from the average predicted value. For example, $Z = 0$ corresponds exactly to the average predicted value and $Z = -2$ means that the measured value is two standard deviations below the average predicted value.

A specific percentile can always be assigned to each Z-score. A Z-score of -1.645 corresponds to the 5% percentile (LLN). If the Z-score is greater than or equal to -1.645 , the measured value is not in the pathological range. The Z-score for the corresponding measured value is marked with a green square in the table of measured values. If the Z-score is smaller than -1.645 , the value is marked with an orange-coloured square (see the guideline on spirometry, Pulmonology. 2015; 69: 146-163).


3) The Miller Quadrant indicates the probability of the existence of a disease and its severity. The relationship between FEV1%VC and FVC is determined and the result is entered in the coordinate system. The coordinate system is divided into the four areas Obstruction, Obstruction & Restriction, Restriction and Normal. The FEV1%VC value is entered as a percentage on the x-axis, the FVC value achieved in comparison to the predicted value is entered as a percentage on the y-axis. The intersection of these values is marked with a cross. The marking crosses are in the colour of the measurement type.

4) The spirometric lung age is determined using the FEV1 value, if FEV1 is calculated depending on age (not for all predicted values). The spirometric lung age is calculated based on the deviation from the predicted value. The spirometric lung age can be displayed via Options, Lung age. On the Examination, Spirometry, Settings, Diagnostic, Parameter screen you can define which predicted value should be used to calculate the spirometric lung age.

5) The following types of unconfirmed reports can be selected: Standard (70% rule for FEV1/FVC and 80% rule for IVC and FVC), COPD-GOLD, clinical or occupational evaluation according to GLI.

6) Table with measured values and predicted values, limit values and the Z-score. Explanation of the evaluation criteria underlying the evaluation functions.

3.5.6 Provocation evaluation




- a** Graphical representation of all the measurements of the measurement series with the selected measured value displayed as a bar, here IVC
- b** Measured value buttons for displaying another measured value in the graphic **a**
- c** Table of measured values – the value selected in the range **b** is highlighted in red
- d** Buttons for selecting a measurement in the measurement series, e.g. to open the selected measurement in the single view (Single Test **g** button)
- e** Environmental data input before the measurement
- f** Comparison between the patient's current and previous evaluation
- g** The selected measurement **d** is displayed as a single test (Same interface as in [3.5.5 Reference and spasmolysis evaluation](#))
- h** Options menu: Print menu, Trend (all evaluations for the patient as a bar diagram and Z-score trend of all evaluations), Export evaluation (e.g. Excel, PDF, e-mail), Change predicted value, Medication Overview, Lung age ¹⁾, Provocation test ²⁾, Reallocate evaluation
- i** Print with system settings
- j** End evaluation

1) The spirometric lung age is determined using the FEV1 value, if FEV1 is calculated depending on age (not for all predicted values). The spirometric lung age is calculated based on the deviation from the predicted value. The spirometric lung age can be displayed via Options, Lung age. On the Examination, Spirometry, Settings, Diagnostic, Parameter screen you can set which predicted value should be used to calculate the spirometric lung age.

2) PD20 provocation dosage: Medication dosage for the 20 percent drop of FEV1 in a provocation measurement compared to the initial value in the reference measurement.

3.5.7 Additional screens of a spirometry evaluation

- **Process control**
(only for reference and spasmolysis evaluations)



- a Flow-volume charts (reference measurement: blue, spasmolysis: orange)
- b Predicted value, in this case GLI
- c Miniature views of the existing measurements
clicking on a measurement displays the respective measurement curve in the a area
- d Table of measured values with predicted values, measured values, Z-score²⁾ and deviations in percentage
- e The selected measurement c is displayed as a single test
- f Options menu: Print menu, Export evaluation (e.g. Excel, PDF, e-mail), Change predicted value, Medication Overview, Reallocate evaluation, repeatability²⁾
- g Print with system settings
- h End evaluation

1) The Z-score indicates by how many standard deviations a certain measured value deviates from the average predicted value. For example, $Z = 0$ corresponds exactly to the average predicted value and $Z = -2$ means that the measured value is two standard deviations below the average predicted value.

A specific percentile can always be assigned to each Z-score. A Z-score of -1.645 corresponds to the 5% percentile (LLN). If the Z-score is greater than or equal to -1.645 , the measured value is not in the pathological range. The Z-score for the corresponding measured value is marked with a green square in the table of measured values. If the Z-score is smaller than -1.645 , the value is marked with an orange-coloured square (see the guideline on spirometry. Pulmonology. 2015; 69: 146-163).


2) In order to make a statement on the quality and plausibility of a measurement series, the FEV1 values of a measurement series and the FVC values of a measurement series are compared with each other. If the deviation is less than 5%, the reproducibility criteria are met. The measured values and limits on which the check is based can be adjusted as required on the Examination, Spirometry, Settings, Menu/Functions screen. Save your input.

This screen is opened via the **Progress Control** button. The **Single Test** button opens the overview again.

The **Process Control** option is used to compare a series of reference and/or spasmolysis measurements in order to check the quality of the patient's cooperation as well as the plausibility of the results. It is a precondition that the measurements have been made in direct succession.

By clicking on the **miniature views of the curves** the measurements can be superimposed. By clicking on the miniature again, the superimposition is disabled again. Strong deviations between the measurement curves show that the patient has not cooperated correctly.



➤ Comparison



With the **Compare** function (at the bottom left of the evaluation overview), the open evaluation can be compared with another evaluation of the patient. The comparison can also be called in the spirometry main menu with **Show Comparison**.

When you click on **Overlap** **a** the measurement curves of both evaluations are superimposed. The **Evaluation** **b** button leads back to the single view of the selected evaluation (pressed date line **c**). The comparative measurement can also be displayed as an individual measurement. For this, click on the date line above the measurement curve **d** and on **Evaluation** **b**. The **arrow buttons** **e** provide you with lists of all the available evaluations for the patient. To open one of these evaluations, select the desired evaluation and click on the **Evaluation** **b** button.

➤ Trend



The trend view is opened in the spirometry main menu with **Show Trend** or in the open evaluation via **Options, Trend**.

The trend view serves to display developments over a longer period of time. All the measurements for a patient are shown as a bar diagram (a) with a table of measured values (b). The selected measured value (c) is displayed in the chart (a) (measured value: blue, predicted value: grey/green). The line with the corresponding measured values is highlighted in colour in the table of measured values (b). To display another measured value, click on the desired measured value button (d).

For measurement series with predicted value GLI: The green areas on the grey bars show the standard value range (e) (upper limit: predicted value, lower limit: LLN). Blue measured value bars whose upper edges are level with the green areas or above can be considered as acceptable. Values below the green areas are in all probability to be considered as pathological. Measured values that are not taken into account with GLI do not have the green standard value ranges for the predicted value display.

In addition, the Z-score values for the patient are shown in the table of measured values (f). Values with a green marking are ≥ -1.645 and therefore acceptable; values with an orange-coloured marking are < -1.645 and probably to be assessed as pathological. The development of the Z-score for the values FEV1, FVC and FEV1/FVC can be displayed via the **ZScore Trend** (g) button³⁾.

3) The Z-score values of the respective measured value are plotted over time (y-axis: Z-score, x-axis: date). The standard value range is within the green lines. Values below the lower green line are to be considered as pathological.

Colour marking of the measured values in the Z-score trend:
 FEV1 — orange
 FVC — pink
 FEV1/FVC — blue

Additional measurements can be displayed using the scroll bar at the bottom of the screen (h). To open a measurement from the trend, select the desired measurement by clicking on the button with the date of creation (i)⁴⁾. Then click on **Evaluation** (i).

4) The trend view also considers the results of already deleted and archived measurements. These results can no longer be shown in the single view.

3.5.8 Confirming the evaluation

Unconfirmed report and report

The unconfirmed report is opened by right-clicking on the evaluation interface. Select **Report** via the context menu. Enter your data in the text field ①. If the **Unconfirmed report or Interpretation** option is selected in the system settings, an automatic system unconfirmed report is already present in the text field¹⁾.

If necessary, older reports can be displayed via the report history (collapsible list above the text input field). When you click on **Confirm** ② your input is saved and the unconfirmed report becomes a (preliminary) report, depending on the report rights of the current user. If your (unconfirmed) report is not yet complete but you want to save it nevertheless without reaching the "Evaluation (pre)confirmed" status, the

report status is reset upon ending (**End**) the evaluation.

Text modules – an aid for writing reports



On the **Examination, Spirometry, Settings, Diagnostic, Report** screen you can create text modules for confirming an evaluation ③. A total of four groups ④ can be stored with up to eight text modules ⑤. The text modules are called in the unconfirmed report dialogue using the keyboard (F5 to F12) ⑥.

A text module can be created from normal text as well as variables. When you use a text module in the unconfirmed report, the actual value from the evaluation is inserted in the report text instead of a variable. The structure of a variable is {VARIABLE} (e.g. FVC predicted value: {FVC_S}). The **Short-cuts for export values** ⑦ button provides you with a list containing all the variables. If the text modules should be shown in the unconfirmed report, make sure that the **Enabled** ⑧ option is enabled. Alternatively, the text modules can be shown in the unconfirmed report via **Options, Texts on** ⑨.

You also have the option of entering a text, which will be automatically shown in each unconfirmed report ⑩. The text can be changed later in the unconfirmed report dialogue. **Save** your input.

1) In addition, further automatic reports may already be included, e.g. according to the criteria clinical, occupational or COPD-GOLD. These automatic reports can be supplemented via **Options, Autom. Report**.

In the settings, you can define which automatic reports should be included in the unconfirmed report by default. To do so, open the **Examination, Spirometry, Settings, Diagnostic, Autom. Report** screen. **Save** your input.



3.5.9 Ending the evaluation

Click on **End** (bottom right) in the evaluation. The End dialogue opens. This is where the "Status of evaluation" ① is defined ¹⁾.

- **evaluation preconfirmed** ②: active if a user with the reporting right "Preconfirm evaluations" has confirmed the unconfirmed report of an evaluation.
- **confirmed** ③ – active if a user with the reporting right "Confirm evaluations" has confirmed the unconfirmed report. The "confirmed" status can be reset if required.
- **printed** ④ – indicates if the evaluation has been printed.
- **indelible** ⑤ – can be selected after reporting has been completed. The evaluation can now only be viewed and can no longer be changed.

Click on **Confirm** ⑥ to close the evaluation.

¹⁾ The assignment of properties (status of evaluation) in the End dialogue makes it easier to find evaluations in the evaluation search.


3.5.10 Optional: Reporting with approval process

If custo diagnostic is used with **approval process**, then authorised persons with the corresponding user rights can save pre-reports of other persons as a report, without having to close the evaluation which was opened previously (shortened workflow) or enter pre-reports/reports directly if the evaluation was created by a person without reporting rights.

The **approval process** is visible in the unconfirmed report dialogue ① of an evaluation. The user or reporter can be changed there (**User name** ②, **Password** ③, **Enter**). During the log-on process, the user rights of the respective user are checked and the software interface is adapted accordingly ④. The reporting is documented in the **evaluation information** ⑤ (context menu).

The **approval process** must be activated user and project-related in the settings and in the custo service centre. The user rights must be set according to the workflow ²⁾. Contact your authorised custo med dealer or custo med.

²⁾ Note: Pre-reporting physicians must have the user right **Preconfirm evaluations**, reporting physicians must have the user rights **Confirm evaluations** and **Change reports of other users**.



3.6 Settings for the spirometry

► Configuring the printout

On the **Examination**, **Spirometry**, **Settings**, **Print**, **General** **a** screen you can define which printed pages are printed when the **Print** button is pressed. In the "Print sequence..." area, select either:

- **Current page** **b** (creates a printout of the open evaluation screen),
- **Pre-set pages** **c** (refer to paragraph below for information on configuration)
- or **Total Printout PA** **d** (professional association, the contents are predefined).

On this screen you can also select the printer for spirometry printouts (on the right half of the screen). **Save** **e** your input.



Defining the contents for the **Pre-set pages** option:

The settings for the printed pages can be found on the **Examination**, **Spirometry**, **Settings**, **Print**, **Printed pages** **f** screen. The contents of the printout can be compiled in the "Printed pages" **g** area.

The **Summary** **h** option is the default setting. This contains a table of measured values, a flow-volume chart, a volume-time curve and an unconfirmed report. The summary can be combined with all other options (report assessment, total table of measured values etc.).

The steps shown are for defining your standard print settings **i**. The standard print settings **i** apply automatically to all other types of printout (print job via the Job Manager, PDF export).

To change the print settings for other types of printout, select the desired printout type in the "Type of printout" area and define the screen contents as described above. **Save** **i** your input.



► Parameters for the spirometry measurement

On the Examination, Spirometry, Settings, Diagnostic, Parameter **a** screen you can set various parameters for the measurement:

Predicted values and area of validity Define which predicted value should be proposed by default for children and adults **b**. If you enable the **Also outside** **c** option in the "Validity" area, the predicted value **b** will also be proposed if the patient data do not match the area of validity of the predicted value. In this case, the predicted values are displayed in brackets in the software interface.

Comparison **d**: You can specify here whether the percentage deviation of the measured values in comparison to the predicted value is to be specified for reference measurements and whether the percentage deviation of the measured values in comparison to the reference or predicted values is to be specified for follow-up measurements.

Spirometric lung age **e**: Select the predicted value according to which the spirometric lung age is to be calculated. The spirometric lung age is determined using the FEV1 value, if FEV1 is calculated depending on age.

Breathing manoeuvre **f**: Adjust the required resting breaths before the breathing manoeuvre.

Identification best value **g**: Define which measured value is to be used to determine the best measurement of a test series. In the "Identification best value" area, select either **Sum of FVC and FEV1**, **FEV1**, **FVC** or **IVC**.

Flow-volume chart **h**: For an ATS-compliant display of the flow-volume chart in custo diagnostic, enable **Show flow-volume chart by ATS...** .


Save your input.

► Display of measured values in the software interface and in the printout

On the **Examination, Spirometry, Settings, Diagnostic, Measur. Display** **a** screen you can set for each predicted value which measured values are to be displayed in the software interface and in the printout (if a different display from the default setting is desired).

In the "Predicted values" area, select the predicted value **b**. Up to seven measured values can then be selected for display **c**. The selected measured values are displayed on the right half of the screen in the "Measurement display" area **d**. The order of the selected measured values can be changed using the **arrow keys** **e**. The measured values IVC, FVC and FEV1 are always displayed and cannot be changed **f**. The measured values IVC, FVC and FEV1 are always displayed and cannot be changed **f**.

By pressing the **reset default values** **g** button, the default settings are displayed again. **Save** **h** your input.



➤ Procedures and functions for the spirometry measurement

These settings can be found on the **Examination**, **Spirometry**, **Settings**, **Menu/Functions**, **Menu/Functions** **a** screen.

Measurement units for the environmental data **b**

You can change the environmental data units in the Environment Data area. The default settings are **°Celsius (°C)**, **hPascal (hPa)** and **Metre (m)**.


Flow control, **Expiratory only** **c**

If this option is enabled, the patient has to exhale forcefully into the device only once for the spirometry measurement (no resting breathing).

Flow control, **Check measurements for repeatability** **d**

In order to make a statement on the quality and plausibility of a measurement series, the FEV1 values of a measurement series and the FVC values of a measurement series are compared with each other. If the deviation is less than 5%, the reproducibility criteria are met. If there are three measurements within the specified limits, a message appears indicating that the measurement series can be ended.

Save **e** your input.




► Quality management

These settings can be found on the Examination, Spirometry, Settings, Menu/Functions, Quality management **a** screen. The Quality management **b** function can be switched on and off as required.

This function supports the correct use of the custo spiro mobile testing device as well as its maintenance and care in order to permanently ensure the quality of the measurements. The checking mechanism examines the best reference measurements of five consecutive patients. If there are deviations from the set limits in five consecutive patients, the system displays corresponding information. The limits are preset as follows:

- If FEV1 is 10 % smaller than the predicted value and the tilt of the measuring head is simultaneously in the red range ($> 10^\circ$) **c**, custo diagnostic displays the message "... Please make sure that the patient is in an upright position during the measurement. ..."
- If IVC or FVC are 10 % greater than the predicted value **d**, custo diagnostic displays the message "... The value for FVC or IVC was more than 10% above the predicted value. We therefore recommend to clean the measuring head and to check the calibration.". The values for checking can be adjusted if necessary.

Save **e** your input.



3.7 Error messages and remedies

Error message: Spirometry device not ready for use

- **Confirm** the error message, close the spirometry software if necessary.
- Disconnect the USB plug of the spirometry testing device from the PC.
- After a few seconds reconnect the device to the PC.
When the LED in the handheld is illuminated, the device is ready for use.
- Call the spirometry software again. This re-initialises the device.
- Afterwards you can continue with the examination.
- If the system still does not recognise your device, custo diagnostic must be restarted.

Error message: Breathing manoeuvre could not be recognised

- If the patient's breathing is too weak or incorrect during the measurement, custo diagnostic might not be able to recognise the breathing manoeuvre.
Without a breathing manoeuvre the system cannot create an evaluation.
- **Confirm** the error message.
- Repeat the measurement, give the patient clear breathing instructions and follow the system's instructions on breathing.
- **Important: The patient must not breathe into the device until you have clicked on the Start button!**

3.8 Appendix

3.8.1 Abbreviations of the spirometry measured values

Abbreviation	Unit	Description
AFEV	l/s	Area under flow volume curve
BF	l/min	Breathing frequency at rest
ERC	l	Expiratory reserve volume
FEF25%FVC	l/s	= MEF75%FVC
FEF25-75%	l/s	Average forced expiratory flow between 25% and 75% of FVC
FEF50%FVC	l/s	= MEF50%FVC
FEF75%FVC	l/s	= MEF25%FVC
FEF75-85%	l/s	Average forced expiratory flow between 75% and 85% of FVC
FET100	s	Forced expiratory time for the total FVC
FET25-75	s	Forced expiratory time between 25 and 75% of FVC
FET95%FVC	s	Forced expiratory time at 95% of FVC
FEV0.5	l	0.5 second capacity (forced expiratory volume in 0.5 seconds)
FEV0.5%FVC	%	Relative 0.5 second capacity of forced vital capacity in percent
FEV0.5/FVC	---	Relative 0.5 second capacity of forced vital capacity
FEV0.5/VC	%	Relative 0.5 second capacity of VC
FEV0.75	l	0.75 second capacity (forced expiratory volume in 0.75 seconds)
FEV0.75/FVC	---	Relative 0.75 second capacity of forced vital capacity
FEV0.75/VC	%	Relative 0.75 second capacity of VC
FEV1	l	Absolute second capacity
FEV1%VC	%	Relative second capacity of vital capacity in percent
FEV1.5	l	1.5 second capacity (forced expiratory volume in 1.5 seconds)
FEV1.5/FVC	---	Relative 1.5 second capacity of forced vital capacity
FEV1.5/VC	%	Relative 1.5 second capacity of VC
FEV1/FEV6	---	Ratio of forced expiratory volume in the first to sixth second
FEV1/FVC	---	Relative second capacity of forced vital capacity
FEV1/VC	%	Relative second capacity of VC
FEV2	l	2 second capacity (forced expiratory volume in 2 seconds)
FEV2/FVC	---	Relative 2 second capacity of forced vital capacity
FEV2/VC	%	Relative 2 second capacity of VC
FEV3	l	3 second capacity (forced expiratory volume in 3 seconds)
FEV3/FVC	---	Relative 3 second capacity of forced vital capacity
FEV3/VC	%	Relative 3 second capacity of VC
FEV6	l	6 second capacity (forced expiratory volume in 6 seconds)
FEV6/FVC	---	Relative 6 second capacity of forced vital capacity
FEV6/VC	%	Relative 6 second capacity of VC
FIF25-75%	l/s	Average forced expiratory flow between 25% and 75% of FVC
FIT100	s	Forced inspiratory time for the total FVC

Abbreviation	Unit	Description
FIV0.5	l	0.5 second capacity (forced inspiratory volume in 0.5 seconds)
FIV0.5/VC	%	Relative 0.5 second capacity of VC (inspirat.) in percent
FIV0.75	l	0.75 second capacity (forced inspiratory volume in 0.75 seconds)
FIV0.75/VC	%	Relative 0.75 second capacity of VC (inspirat.) in percent
FIV1	l	Second capacity (forced inspiratory volume in 1 second)
FIV1%VC	%	Relative second capacity of VC (inspirat.) in percent
FIV1.5	l	1.5 second capacity (forced inspiratory volume in 1.5 seconds)
FIV1.5/VC	%	Relative 1.5 second capacity of VC (inspirat.) in percent
FIV2	l	2 second capacity (forced inspiratory volume in 2 seconds)
FIV2/VC	%	Relative 2 second capacity of VC (inspirat.) in percent
FIV3	l	3 second capacity (forced inspiratory volume in 3 seconds)
FIV3/VC	%	Relative 3 second capacity of VC (inspirat.) in percent
FVC	l	Forced vital capacity
IC	[l]	Inspiratory capacity (amount of air that can be inhaled after normal expiration) IRC + T
IRC	l	Inspiratory reserve volume
IVC	l	Inspiratory vital capacity
Lung age	Years	The patient's spirometric lung age is determined using the measured FEV1, if FEV1 is calculated depending on age (different depending on the predicted value).
MEF25%FVC	l/s	Forced expiratory flow at 25% of FVC
MEF50%FVC	l/s	Forced expiratory flow at 50% of FVC
MEF75%FVC	l/s	Forced expiratory flow at 75% of FVC
MIF25%FVC	l/s	Forced inspiratory flow at 25% of FVC
MIF50%FVC	l/s	Forced inspiratory flow at 50% of FVC
MIF75%FVC	l/s	Forced inspiratory flow at 75% of FVC
MVV	l	Maximum voluntary ventilation
OBQ		Obstruction ratio
PEF	l/s	Peak flow, maximum expiratory flow
PIF	l/s	Peak flow, maximum inspiratory flow
tE	s	Average time of expiration at rest
tI	s	Average time of inspiration at rest
TV	l	Tidal Volume
VCmax	l	Maximum vital capacity, inspiratory or expiratory
VTtI	l/s	Average inspiratory flow at rest

3.8.2 Calculation tables for predicted values

Multicèntrico di Barcelona	39	ECCS/Quanjer	43
Polgar79	40	Hankinson	44
Polgar71	40	HSU	48
Crapo	40	Schindl	48
Morris	41	ECCS-children/Quanjer	48
Austrian reference values	41	Zapletal	49
Cherniak	41	GLI - Global Lung Function Initiative	49
Knudson	42	Kainu, Finland	50
Ulmer	43	Siriraj, Thailand	50
Baur	43	Danish Reference Values	50

The predicted values define their areas of validity using age, height, weight, ethnicity etc. A suitable predicted value for the measurement is allocated to the patient according to his/her data. The standard setting for children and adults is GLI.

Abbreviations in the calculation tables

- A = Age
- H = Height
- G = Weight
- B = Broca index = Weight: (Height - 100)
- Fi = Obesity = $H : \sqrt[3]{W}$
- M = Predicted value
- S = Coefficient of variation
- AfrAm = Afro-American ethnicity
- NEAsia = North-East Asian ethnicity
- SEAsia = South-East Asian ethnicity
- Other = other ethnic groups

Multicèntrico di Barcelona		Boys	Girls
		H = 85 – 180 cm A = 6 – 20 years	H = 85 – 180 cm A = 6 – 20 years
FVC	[l]	$0.02800 * H + 0.03451 * G + 0.05728 * A - 3.21$	$0.03049 * H + 0.02220 * G + 0.03550 * A - 3.04$
FEV1	[l]	$0.02483 * H + 0.02266 * G + 0.07148 * A - 2.91$	$0.02866 * H + 0.01713 * G + 0.02955 * A - 2.87$
MVV	[l]	$(0.02483 * H + 0.02266 * G + 0.07148 * A - 2.91) * 37.5$	$(0.02866 * H + 0.01713 * G + 0.02955 * A - 2.87) * 37.5$
FEF25-75%	[l/s]	$0.038 * H + 0.140 * A - 4.33$	$0.046 * H + 0.051 * A - 4.30$
PEF	[l/s]	$0.075 * H + 0.275 * A - 9.08$	$0.073 * H + 0.134 * A - 7.57$
MEF25%FVC	[l/s]	$0.024 * H + 0.066 * A - 2.61$	$0.027 * H + 0.032 * A - 2.68$
MEF50%FVC	[l/s]	$0.017 * H + 0.157 * A + 0.029 * G - 2.17$	$0.046 * H + 0.067 * A - 4.17$
FEV1%VC	[%]	$-0.1902 * A + 85.58$	$-0.224 * A - 0.1126 * G + 94.88$
		Men	Women
		H = 150 – 200 cm A = > 20 years	H = 150 – 200 cm A = > 20 years
FVC	[l]	$0.0678 * H - 0.0147 * A - 6.05$	$0.0454 * H - 0.0221 * A - 2.83$
FEV1	[l]	$0.0499 * H - 0.0211 * A - 3.84$	$0.0317 * H - 0.0250 * A - 1.23$
MVV	[l]	$(0.0499 * H - 0.0211 * A - 3.84) * 37.5$	$(0.0317 * H - 0.0250 * A - 1.23) * 37.5$
FEF25-75%	[l/s]	$0.0392 * H - 0.0430 * A - 1.16$	$0.0230 * H - 0.0456 * A - 1.11$
PEF	[l/s]	$0.0945 * H - 0.0209 * A - 5.77$	$0.0448 * H - 0.0304 * A - 0.35$
MEF25%FVC	[l/s]	$0.0190 * H - 0.0356 * A - 0.14$	$0.02 * H - 0.031 * A - 0.0062 * G - 0.21$
MEF50%FVC	[l/s]	$0.0517 * H - 0.0397 * A - 2.40$	$0.0242 * H - 0.0418 * A - 1.62$
FEV1%VC	[%]	$-0.1902 * A + 85.58$	$-0.224 * A - 0.1126 * G + 94.88$

Polgar79		Boys H = 85 – 180 cm A = 4 – 17 years	Girls H = 85 – 180 cm A = 4 – 18 years
		Men H = 150 – 200 cm A = 18 – 120 years	Women H = 150 – 200 cm A = 18 – 120 years
FVC	[l]	$2.12 * 0.000001 * H^{2.81}$	$2.34 * 0.000001 * H^{2.78}$
IVC	[l]	$2.12 * 0.000001 * H^{2.81}$	$2.34 * 0.000001 * H^{2.78}$
FEF25-75%	[l/s]	$(219.66 + 2.72 * H) : 60$	$(219.66 + 2.72 * H) : 60$
PEF	[l/s]	$(467.96 + 5.59 * H) : 60$	$(376.51 + 4.85 * H) : 60$

Polgar71		Boys H = 85 – 180 cm A = 4 – 17 years	Girls H = 85 – 180 cm A = 4 – 18 years
		Men H = 150 – 200 cm A = 18 – 120 years	Women H = 150 – 200 cm A = 18 – 120 years
FVC	[l]	$4.4 * 0.000001 * H^{2.67}$	$3.3 * 0.000001 * H^{2.72}$
IVC	[l]	$4.4 * 0.000001 * H^{2.67}$	$3.3 * 0.000001 * H^{2.72}$
FEV1	[l]	$2.1 * 0.000001 * H^{2.8}$	$2.1 * 0.000001 * H^{2.8}$
MVV	[l]	$99.507 + 1.276 * H$	$99.507 + 1.276 * H$
FEF25-75%	[l/s]	$(207.7 + 2.621 * H) : 60$	$(207.7 + 2.621 * H) : 60$
PEF	[l/s]	$(425.5714 + 5.2428 * H) : 60$	$(-425.5714 + 5.2428 * H) : 60$
PIF	[l/s]	$5.26 + 0.06 * H$	$5.26 + 0.06 * H$

Crapo		Men H = 150 – 220 cm A = 18 – 120 years	Women H = 150 – 220 cm A = 18 – 120 years
FVC	[l]	$6.00 * H - 0.0214 * A - 4.650$	$4.91 * H - 0.0216 * A - 3.590$
IVC	[l]	$6.00 * H - 0.0214 * A - 4.650$	$4.91 * H - 0.0216 * A - 3.590$
FEV0.5	[l]	$3.27 * H - 0.0152 * A - 1.914$	$2.38 * H - 0.0185 * A - 0.809$
FEV1	[l]	$4.14 * H - 0.0244 * A - 2.190$	$3.42 * H - 0.0255 * A - 1.578$
MVV	[l]	$(4.14 * H - 0.0244 * A - 2.190) * 37.5$	$(3.42 * H - 0.0255 * A - 1.578) * 37.5$
FEV3	[l]	$5.35 * H - 0.0271 * A - 3.512$	$4.42 * H - 0.0257 * A - 2.745$
FEV1%VC	[%]	$13.0 * H - 0.152 * A + 110.49$	$20.20 * H - 0.252 * A + 126.58$
FEV3/VC	[%]	$6.27 * H - 0.145 * A + 112.09$	$9.37 * H - 0.163 * A + 118.16$
FEF25-75%	[l/s]	$2.04 * H - 0.038 * A + 2.133$	$1.54 * H - 0.046 * A + 2.683$

Morris		Men H = 150 – 220 cm A = 20 – 120 years	Women H = 150 – 220 cm A = 20 – 120 years
FVC	[l]	$5.83 * H - 0.025 * A - 4.241$	$4.52 * H - 0.024 * A - 2.852$
IVC	[l]	$5.83 * H - 0.025 * A - 4.241$	$4.52 * H - 0.024 * A - 2.852$
FEV1	[l]	$3.62 * H - 0.032 * A - 1.260$	$3.50 * H - 0.025 * A - 1.932$
MVV	[l]	$3.62 * H - 0.032 * A - 1.260 * 37.5$	$3.50 * H - 0.025 * A - 1.932 * 37.5$
FEV1%VC	[%]	$107.12 - 12.28 * H - 0.2422 * A$	$88.70 - 2.67 * H - 0.1815 * A$
FEF25–75%	[l/s]	$1.85 * H - 0.045 * A + 2.513$	$2.36 * H - 0.030 * A + 0.551$

Austrian reference values ¹⁾		Men H = 1.44 – 2.00 m A = 18 – 90 years	Women H = 1.40 – 1.90 m A = 16 – 90 years
FVC	[l]	$-11.606 + 8.172H - 0.0339A * H + 1.2869 \ln(A)$	$-10.815 + 6.640H - 0.0408A * H + 1.7293 \ln(A)$
FEV1	[l]	$-8.125 + 6.212H - 0.0300A * H + 0.9770 \ln(A)$	$-6.995 + 5.174 - 0.0314A * H + 1.0251 \ln(A)$
PEF	[l/s]	$(1.798 + 2.311 \ln(H) + 0.0159A - 0.000248A^2)^2$	$(1.832 + 1.838 \ln(H) + 0.0078A - 0.000172A^2)^2$
MEF75%FVC	[l/s]	$(1.581 + 1.854 \ln(H) + 0.0213A - 0.000283A^2)^2$	$(1.779 + 1.421 \ln(H) + 0.0096A - 0.000179A^2)^2$
MEF50%FVC	[l/s]	$(1.490 + 1.290 \ln(H) + 0.0125A - 0.000218A^2)^2$	$(1.561 + 1.177 \ln(H) + 0.0045A - 0.000140A^2)^2$
MEF25%FVC	[l/s]	$(1.314 + 0.898 \ln(H) - 0.0083A - 0.000026A^2)^2$	$(1.372 + 0.938 \ln(H) - 0.0152A + 0.000036A^2)^2$
FEV1%VC	[%]	$101.99 - 1.191H^2 - 3.962 \ln(A)$	$118.993 - 3.0320H^2 - 6.9053 \ln(A)$
		Boys H = 1.09 – 1.96 m A = 5 – 17.99 years	Girls H = 1.10 – 1.82 m A = 5 – 15.99 years
FVC	[l]	$\exp(-1.142 + 1.259H + 0.004070A \sqrt{W})$	$\exp(-3.842 + 4.1632 \sqrt{H} + 0.1341 \sqrt{A} - 1.614Fi)$
FEV1	[l]	$\exp(-1.178 + 1.221H + 0.003841A \sqrt{W})$	$\exp(-3.877 + 3.9808 \sqrt{H} + 0.1485 \sqrt{A} - 1.322Fi)$
PEF	[l/s]	$\exp(-0.214 + 0.921H + 0.0467A + 0.0020W)$	$\exp(0.411 + 1.793 \ln(H) + 0.4251 \ln(A) - 0.910Fi)$
MEF75%FVC	[l/s]	$\exp(-0.077 + 0.770H + 0.0373A + 0.0025W)$	$\exp(0.455 + 1.616 \ln(H) + 0.3738 \ln(A) - 0.861Fi)$
MEF50%FVC	[l/s]	$\exp(-0.522 + 0.843H + 0.0300A + 0.0035W)$	$\exp(0.256 + 1.643 \ln(H) + 0.3481 \ln(A) - 1.089Fi)$
MEF25%FVC	[l/s]	$\exp(-1.576 + 1.166H + 0.0219A + 0.0021W)$	$\exp(-0.772 + 2.002 \ln(H) + 0.3063 \ln(A) - 0.409Fi)$
FEV1%VC	[%]	$(101.99 - 1.191H^2 - 3.962 \ln(A))$	92

Cherniak ²⁾		Men H = 150 – 190 cm A = 15 – 79 years	Women H = 150 – 190 cm A = 15 – 79 years
FVC	[ml]	$47.6 * H - 14 * A - 3180$	$30.7 * H - 15 * A - 1310$
FEV1	[ml]	$35.9 * H - 23 * A - 1510$	$23.7 * H - 19 * A - 0190$
MVV	[ml]	$(35.9 * H - 23 * A - 1510) * 37.5$	$(23.7 * H - 19 * A - 0190) * 37.5$
PEF	[ml/s]	$57.6 * H - 24 * A + 0230$	$35.9 * H - 18 * A + 1130$
MEF75%FVC	[ml/s]	$35.6 * H - 20 * A + 2730$	$27.1 * H - 19 * A + 2150$
MEF50%FVC	[ml/s]	$25.7 * H - 30 * A + 2400$	$24.5 * H - 23 * A + 1430$
MEF25%FVC	[ml/s]	$14.1 * H - 41 * A + 1610$	$09.2 * H - 35 * A + 2220$
		Boys H = 75 – 180 cm A = 3 – 17 years	Girls H = 75 – 180 cm A = 3 – 17 years
FVC	[ml]	$40.53 * H + 51.34 * A - 3655$	$27.86 * H + 90.96 * A - 2554$

Knudson ³⁾		Men H = 150 – 195 cm A = 25 – 80 years	Women H = 150 – 195 cm A = 25 – 80 years
FVC	[ml]	$65 * H - 29 * A - 5460$	$37 * H - 22 * A - 1770$
FEV1	[ml]	$52 * H - 27 * A - 4200$	$27 * H - 21 * A - 790$
MVV	[ml]	$(52 * H - 27 * A - 4200) * 37.5$	$(27 * H - 21 * A - 790) * 37.5$
FEV1%VC	[%]	$0.087 * H - 0.14 * A + 103.64$	$0.111 * H - 0.109 * A + 107.38$
PEF	[ml/s]	$94 * H - 35 * A - 5993$	$49 * H - 25 * A - 735$
MEF75%FVC	[ml/s]	$88 * H - 35 * A - 5620$	$43 * H - 25 * A - 130$
MEF50%FVC	[ml/s]	$69 * H - 15 * A - 5400$	$35 * H - 13 * A - 440$
MEF25%FVC	[ml/s]	$44 * H - 12 * A - 4140$	$-14 * A + 3040$
		Boys H = 140 – 193 cm A = 12 – 25 years	Girls H = 140 – 193 cm A = 12 – 25 years
FVC	[ml]	$59.0 * H - 73.9 * A - 6887$	$41.6 * H + 69.9 * A - 4447$
FEV1	[ml]	$51.9 * H - 6118$	$35.1 * H + 6.94 * A - 3762$
MVV	[ml]	$(51.9 * H - 6118) * 37.5$	$(35.1 * H + 6.94 * A - 3762) * 37.5$
FEV1%VC	[%]	$-0.0813 * H + 100.439$	$-0.1909 * H + 0.6655 * A + 109.97$
PEF	[ml/s]	$78.0 * H + 166 * A - 8060$	$49.0 * H + 157 * A - 3916$
MEF75%FVC	[ml/s]	$70.0 * H + 147 * A - 7054$	$44.0 * H + 144 * A - 3365$
MEF50%FVC	[ml/s]	$54.3 * H + 115 * A - 6385$	$28.8 * H + 111 * A - 2304$
MEF25%FVC	[ml/s]	$39.7 * H - 5.7 * A - 4242$	$24.3 * H + 292.3 * A - 7.5 * A^2 - 4400.9$
		Boys H = 112 – 155 cm A = 6 – 12 years	Girls H = 112 – 155 cm A = 6 – 12 years
FVC	[ml]	$40.9 * H - 3376$	$43.0 * H - 3749$
FEV1	[ml]	$34.0 * H - 2814$	$33.6 * H - 2758$
MVV	[ml]	$(34.0 * H - 2814) * 37.5$	$(33.6 * H - 2758) * 37.5$
FEV1%VC	[%]	$0.0813 * H + 100.439$	$-0.1909 * H + 0.6655 * A + 109.97$
PEF	[ml/s]	$78.0 * H + 166 * A - 8060$	$49.9 * H + 157 * A - 3916$
MEF75%FVC	[ml/s]	$70.0 * H + 147 * A - 7054$	$44.0 * H + 144 * A - 3365$
MEF50%FVC	[ml/s]	$37.8 * H + 2545$	$184.6 * A + 736$
MEF25%FVC	[ml/s]	$17.1 * H - 1014.9$	$10.9 * H - 165.7$

Ulmer ⁴⁾		Men H = 150 – 195 cm A = 15 – 75 years G = 40 – 170 kg	Women H = 150 – 195 cm A = 15 – 75 years G = 40 – 170 kg
IVC	[ml]	$82.243 * H - 20.4 * A - 8420.5 - 69.8 * B$	$56.695 * H - 19.4 * A - 5096 - 69.7 * B$
IRC	[ml]	$47.291 * H - 11.3 * A - 6632 + 1297.3 * B$	$35.751 * H - 6.4 * A - 4241.4 - 1016.1 * B$
ERC	[ml]	$41.995 * H - 7.8 * A - 3523.8 - 1875 * B$	$12.126 * H - 14.4 * A + 136 - 624.6 * B$
FVC	[ml]	$77.576 * H - 21.7 * A - 7769.5 - 151.3 * B$	$52.467 * H - 19.9 * A - 4412.3 - 400.4 * B$
FEV1	[ml]	$53.212 * H - 26.1 * A - 4234 - 71.8 * B$	$23.939 * H - 20.7 * A - 641.6 - 209 * B$
MVV	[ml]	$(53.212 * H - 26.1 * A - 4234 - 71.8 * B) * 37.5$	$(23.939 * H - 20.7 * A - 641.6 - 209 * B) * 37.5$
PEF	[ml/s]	$66.067 * H - 20.8 * A - 2981.3 - 1249.3 * B$	$55.175 * H - 31.4 * A - 1683.4 - 115.1 * B$
MEF50%FVC	[ml/s]	$30.584 * H - 44 * A + 672.3 + 668.5 * B$	$26.181 * H - 22.4 * A + 2618.1 + 124 * B$
MEF25%FVC	[ml/s]	$25.108 * H - 39 * A - 1254.2 + 697.4 * B$	$20.129 * H - 35.2 * A - 438.6 + 593.6 * B$

Baur ⁵⁾		Men H = 1.55 – 1.95 m A = 18 – 70 years	Women H = 1.45 – 1.80 m A = 18 – 70 years
FVC	[l]	$6.00 * H - 0.0214 * A - 4.650$	$4.91 * H - 0.0216 * A - 3.590$
FEV1	[l]	$4.14 * H - 0.0244 * A - 2.190$	$3.42 * H - 0.0255 * A - 1.578$
MVV	[l]	$(4.14 * H - 0.0244 * A - 2.190) * 37.5$	$(3.42 * H - 0.0255 * A - 1.578) * 37.5$
PEF	[l/s]	$6.14 * H - 0.043 * A + 0.15$	$5.50 * H - 0.030 * A - 1.11$
MEF75%FVC	[l/s]	$5.46 * H - 0.029 * A - 0.47$	$3.22 * H - 0.025 * A + 1.60$
MEF50%FVC	[l/s]	$3.79 * H - 0.031 * A - 0.35$	$2.45 * H - 0.025 * A + 1.16$
MEF25%FVC	[l/s]	$2.61 * H - 0.026 * A - 1.34$	$1.05 * H - 0.025 * A + 1.11$

ECCS/Quanjer ⁶⁾		Men H = 150 – 195 cm A = 25 – 75 years	Women H = 150 – 190 cm A = 25 – 75 years
IVC	[ml]	$61.03 * H - 28 * A - 4654$	$46.64 * H - 26 * A - 3284$
FVC	[ml]	$57.57 * H - 26 * A - 4345$	$44.26 * H - 26 * A - 2887$
FEV1	[ml]	$43.01 * H - 29 * A - 2492$	$39.53 * H - 25 * A - 2604$
MVV	[ml]	$(43.01 * H - 29 * A - 2492) * 37.5$	$(39.53 * H - 25 * A - 2604) * 37.5$
FEV1%VC	[%]	$87.21 * H - 0.179 * A$	$89.10 * H - 0.192 * A$
PEF	[ml/s]	$61.46 * H - 43 * A + 154$	$55.01 * H - 30 * A - 1106$
MEF75%FVC	[ml/s]	$54.59 * H - 29 * A - 470$	$32.18 * H - 25 * A + 1596$
MEF50%FVC	[ml/s]	$37.94 * H - 31 * A - 352$	$24.50 * H - 25 * A + 1156$
MEF25%FVC	[ml/s]	$26.05 * H - 26 * A - 1336$	$10.50 * H - 25 * A + 1107$
FEF25–75%	[ml/s]	$19.4 * H - 43.0 * A + 2700.0$	$12.5 * H - 34.0 * A + 2920.0$

Hankinson		Boys
		Caucasian/Asian
		H = 75 – 180 cm A = 4 – 19 years
FEV1%VC	[%]	$88.066 + (-0.2066 * A)$
FEV1	[l]	$0.7453 + (-0.04106 * A) + (0.004477 * A * A) + (0.00014098 * H * H)$
MVV	[l]	$(0.7453 + (-0.04106 * A) + (0.004477 * A * A) + (0.00014098 * H * H)) * 37.5$
FEV6	[l]	$-0.3119 + (-0.18612 * A) + (0.009717 * A * A) + (0.00018188 * H * H)$
FVC	[l]	$-0.2584 + (-0.20415 * A) + (0.010133 * A * A) + (0.00018642 * H * H)$
IVC	[l]	$-0.2584 + (-0.20415 * A) + (0.010133 * A * A) + (0.00018642 * H * H)$
PEF	[l/s]	$-0.5962 + (-0.12357 * A) + (0.013135 * A * A) + (0.00024962 * H * H)$
FEF25–75%	[l/s]	$-1.0863 + (0.13939 * A) + (0.00010345 * H * H)$
		Afro-American
		H = 75 – 180 cm A = 4 – 19 years
FEV1%VC	[%]	$89.239 + (-0.1828 * A)$
FEV1	[l]	$-0.7048 + (-0.05711 * A) + (0.004316 * A * A) + (0.00013194 * H * H)$
MVV	[l]	$((-0.7048 + (-0.05711 * A) + (0.004316 * A * A) + (0.00013194 * H * H)) * 37.5$
FEV6	[l]	$-0.5525 + (-0.14107 * A) + (0.007241 * A * A) + (0.00016429 * H * H)$
FVC	[l]	$-0.4971 + (-0.15497 * A) + (0.007701 * A * A) + (0.00016643 * H * H)$
IVC	[l]	$-0.4971 + (-0.15497 * A) + (0.007701 * A * A) + (0.00016643 * H * H)$
PEF	[l/s]	$-0.2684 + (-0.28016 * A) + (0.018202 * A * A) + (0.00027333 * H * H)$
FEF25–75%	[l/s]	$-1.1627 + (0.12314 * A) + (0.00010461 * H * H)$
		Latin-American
		H = 75 – 180 cm A = 4 – 19 years
FEV1%VC	[%]	$90.024 + (-0.2186 * A)$
FEV1	[l]	$-0.8218 + (-0.04248 * A) + (0.004291 * A * A) + (0.00015104 * H * H)$
MVV	[l]	$(-0.8218 + (-0.04248 * A) + (0.004291 * A * A) + (0.00015104 * H * H)) * 37.5$
FEV6	[l]	$-0.6646 + (-0.11270 * A) + (0.007306 * A * A) + (0.00017840 * H * H)$
FVC	[l]	$-0.7571 + (-0.09520 * A) + (0.006619 * A * A) + (0.00017823 * H * H)$
IVC	[l]	$-0.7571 + (-0.09520 * A) + (0.006619 * A * A) + (0.00017823 * H * H)$
PEF	[l/s]	$-0.9537 + (-0.19602 * A) + (0.014497 * A * A) + (0.00030243 * H * H)$
FEF25–75%	[l/s]	$-1.3592 + (0.10529 * A) + (0.00014473 * H * H)$

Hankinson		Girls
		Caucasian/Asian
		H = 75 – 180 cm A = 4 – 17 years
FEV1%VC	[%]	$90.809 + (-0.2125 * A)$
FEV1	[l]	$-0.8710 + (0.06537 * A) + (0.00011496 * H * H)$
MVV	[l]	$(-0.8710 + (0.06537 * A) + (0.00011496 * H * H)) * 37.5$
FEV6	[l]	$-1.1925 + (0.06544 * A) + (0.00014395 * H * H)$
FVC	[l]	$-1.2082 + (0.05916 * A) + (0.00014815 * H * H)$
IVC	[l]	$-1.2082 + (0.05916 * A) + (0.00014815 * H * H)$
PEF	[l/s]	$-3.6181 + (0.60644 * A) + (-0.016846 * A * A) + (0.00018623 * H * H)$
FEF25–75%	[l/s]	$-2.5284 + (0.52490 * A) + (-0.015309 * A * A) + (0.00006982 * H * H)$
		Afro-American
		H = 75 – 180 cm A = 4 – 17 years
FEV1%VC	[%]	$91.655 + (-0.2039 * A)$
FEV1	[l]	$-0.9630 + (0.05799 * A) + (0.00010846 * H * H)$
MVV	[l]	$((-0.9630 + (0.05799 * A) + (0.00010846 * H * H)) * 37.5$
FEV6	[l]	$0.6370 + (-0.04243 * A) + (0.003508 * A * A) + (0.00013497 * H * H)$
FVC	[l]	$-0.6166 + (-0.04687 * A) + (0.003602 * A * A) + (0.00013606 * H * H)$
IVC	[l]	$-0.6166 + (-0.04687 * A) + (0.003602 * A * A) + (0.00013606 * H * H)$
PEF	[l/s]	$-1.2398 + (0.16375 * A) + (0.00019746 * H * H)$
FEF25–75%	[l/s]	$-2.5379 + (0.43755 * A) + (-0.012154 * A * A) + (0.00008572 * H * H)$
		Latin-American
		H = 75 – 180 cm A = 4 – 17 years
FEV1%VC	[%]	$92.360 + (-0.2248 * A)$
FEV1	[l]	$-0.9641 + (0.06490 * A) + (0.00012154 * H * H)$
MVV	[l]	$((-0.9641 + (0.06490 * A) + (0.00012154 * H * H)) * 37.5$
FEV6	[l]	$-1.2410 + (0.07625 * A) + (0.00014106 * H * H)$
FVC	[l]	$-1.2507 + (0.07501 * A) + (0.00014246 * H * H)$
IVC	[l]	$-1.2507 + (0.07501 * A) + (0.00014246 * H * H)$
PEF	[l/s]	$-3.2549 + (0.47495 * A) + (-0.013193 * A * A) + (0.00022203 * H * H)$
FEF25–75%	[l/s]	$-2.1825 + (0.42451 * A) + (-0.012415 * A * A) + (0.00009610 * H * H)$

Hankinson		Men
		Caucasian/Asian
		H = 150 – 200 cm A = 20 – 120 years
FEV1%VC	[%]	$88.066 + (-0.2066 * A)$
FEV1	[l]	$0.5536 + (-0.01303 * A) + (-0.000172 * A * A) + (0.00014098 * H * H)$
MVV	[l]	$(0.5536 + (-0.01303 * A) + (-0.000172 * A * A) + (0.00014098 * H * H)) * 37.5$
FEV6	[l]	$0.1102 + (-0.00842 * A) + (-0.000223 * A * A) + (0.00018188 * H * H)$
FVC	[l]	$-0.1933 + (0.00064 * A) + (-0.000269 * A * A) + (0.00018642 * H * H)$
IVC	[l]	$-0.1933 + (0.00064 * A) + (-0.000269 * A * A) + (0.00018642 * H * H)$
PEF	[l/s]	$1.0523 + (0.08272 * A) + (-0.001301 * A * A) + (0.00024962 * H * H)$
FEF25–75%	[l/s]	$2.7006 + (-0.04995 * A) + (0.00010345 * H * H)$
		Afro-American
		H = 150 – 200 cm A = 20 – 120 years
FEV1%VC	[%]	$89.239 + (-0.1828 * A)$
FEV1	[l]	$0.3411 + (-0.02309 * A) + (0.00013194 * H * H)$
MVV	[l]	$(0.3411 + (-0.02309 * A) + (0.00013194 * H * H)) * 37.5$
FEV6	[l]	$-0.0547 + (-0.02114 * A) + (0.00016429 * H * H)$
FVC	[l]	$-0.1517 + (-0.01821 * A) + (0.00016643 * H * H)$
IVC	[l]	$-0.1517 + (-0.01821 * A) + (0.00016643 * H * H)$
PEF	[l/s]	$2.2257 + (-0.04082 * A) + (0.00027333 * H * H)$
FEF25–75%	[l/s]	$2.1477 + (-0.04238 * A) + (0.00010461 * H * H)$
		Latin-American
		H = 150 – 200 cm A = 20 – 120 years
FEV1%VC	[%]	$90.024 + (-0.2186 * A)$
FEV1	[l]	$0.6306 + (-0.02928 * A) + (0.00015104 * H * H)$
MVV	[l]	$(0.6306 + (-0.02928 * A) + (0.00015104 * H * H)) * 37.5$
FEV6	[l]	$0.5757 + (-0.02860 * A) + (0.00017840 * H * H)$
FVC	[l]	$0.2376 + (-0.00891 * A) + (-0.000182 * A * A) + (0.00017823 * H * H)$
IVC	[l]	$0.2376 + (-0.00891 * A) + (-0.000182 * A * A) + (0.00017823 * H * H)$
PEF	[l/s]	$0.0870 + (0.06580 * A) + (-0.001195 * A * A) + (0.00030243 * H * H)$
FEF25–75%	[l/s]	$1.7503 + (-0.05018 * A) + (0.00014473 * H * H)$

Hankinson		Women
		Caucasian/Asian
		H = 140 – 200 cm A = 18 – 120 years
FEV1%VC	[%]	$90.809 + (-0.2125 * A)$
FEV1	[l]	$0.4333 + (-0.00361 * A) + (-0.000194 * A * A) + (0.00011496 * H * H)$
MVV	[l]	$(0.4333 + (-0.00361 * A) + (-0.000194 * A * A) + (0.00011496 * H * H)) * 37.5$
FEV6	[l]	$-0.1373 + (0.01317 * A) + (-0.000352 * A * A) + (0.00014395 * H * H)$
FVC	[l]	$-0.3560 + (0.01870 * A) + (-0.000382 * A * A) + (0.00014815 * H * H)$
IVC	[l]	$-0.3560 + (0.01870 * A) + (-0.000382 * A * A) + (0.00014815 * H * H)$
PEF	[l/s]	$0.9267 + (0.06929 * A) + (-0.001031 * A * A) + (0.00018623 * H * H)$
FEF25–75%	[l/s]	$2.3670 + (-0.01904 * A) + (-0.000200 * A * A) + (0.00006982 * H * H)$
		Afro-American
		H = 140 – 200 cm A = 18 – 120 years
FEV1%VC	[%]	$91.655 + (-0.2039 * A)$
FEV1	[l]	$0.3433 + (-0.01283 * A) + (-0.000097 * A * A) + (0.00010846 * H * H)$
MVV	[l]	$(0.3433 + (-0.01283 * A) + (-0.000097 * A * A) + (0.00010846 * H * H)) * 37.5$
FEV6	[l]	$-0.1981 + (0.00047 * A) + (-0.000230 * A * A) + (0.00013497 * H * H)$
FVC	[l]	$-0.3039 + (0.00536 * A) + (-0.000265 * A * A) + (0.00013606 * H * H)$
IVC	[l]	$-0.3039 + (0.00536 * A) + (-0.000265 * A * A) + (0.00013606 * H * H)$
PEF	[l/s]	$1.3597 + (0.03458 * A) + (-0.000847 * A * A) + (0.00019746 * H * H)$
FEF25–75%	[l/s]	$2.0828 + (-0.03793 * A) + (0.00008572 * H * H)$
		Latin-American
		H = 140 – 200 cm A = 18 – 120 years
FEV1%VC	[%]	$92.360 + (-0.2248 * A)$
FEV1	[l]	$0.4529 + (-0.01178 * A) + (-0.000113 * A * A) + (0.00012154 * H * H)$
MVV	[l]	$(0.4529 + (-0.01178 * A) + (-0.000113 * A * A) + (0.00012154 * H * H)) * 37.5$
FEV6	[l]	$0.2033 + (0.00020 * A) + (-0.000232 * A * A) + (0.00014106 * H * H)$
FVC	[l]	$0.1210 + (0.00307 * A) + (-0.000237 * A * A) + (0.00014246 * H * H)$
IVC	[l]	$0.1210 + (0.00307 * A) + (-0.000237 * A * A) + (0.00014246 * H * H)$
PEF	[l/s]	$0.2401 + (0.06174 * A) + (-0.001023 * A * A) + (0.00022203 * H * H)$
FEF25–75%	[l/s]	$1.7456 + (-0.01195 * A) + (-0.000291 * A * A) + (0.00009610 * H * H)$

HSU		Boys H = 75 – 180 cm A = 7 – 18 years	Girls H = 75 – 180 cm A = 7 – 18 years
FVC	[l]	$(3.58 : 10000) * H^{3.18} : 1000$	$(2.57 : 1000) * H^{2.78} : 1000$
IVC	[l]	$(3.58 : 10000) * H^{3.18} : 1000$	$(2.57 : 1000) * H^{2.78} : 1000$
FEV1	[l]	$(7.74 : 10000) * H^{3.00} : 1000$	$(3.79 : 1000) * H^{2.68} : 1000$
MVV	[l]	$(7.74 : 10000) * H^{3.00} : 1000 * 37.5$	$(3.79 : 1000) * H^{2.68} : 1000 * 37.5$
PEF	[l/s]	$((3.35 : 10000) * H^{2.79}) : 60$	$((2.58 : 1000) * H^{2.37}) : 60$
FEF25–75%	[l/s]	$((7.98 : 10000) * H^{2.46}) : 60$	$((3.79 : 1000) * H^{2.16}) : 60$

Schindl ^[7]		Boys H = 110 – 180 cm A = 10 – 18 years	Girls H = 110 – 180 cm A = 10 – 18 years
FVC	[ml]	$49.2 * H + 118.2 * A - 6006.0$	$41.7 * H + 91.3 * A - 4660.6$
FEV1	[ml]	$41.9 * H + 79.0 * A - 4674.4$	$41.9 * H + 70.6 * A - 4176.1$
PEF	[ml/s]	$76.8 * H + 224.2 * A - 8381.5$	$62.1 * H + 176.3 * A - 5623.2$
MEF75%FVC	[ml/s]	$56.2 * H + 175.4 * A - 5530.3$	$46.5 * H + 154.7 * A - 3627.9$
MEF50%FVC	[ml/s]	$41.5 * H + 109.5 * A - 3988.0$	$48.3 * H + 115.6 * A - 4896.6$
MEF25%FVC	[ml/s]	$30.3 * H + 39.0 * A - 3059.9$	$38.8 * H + 51.4 * A - 4331.9$

ECCS–children/Quanjer		Boys H = 75 – 180 cm A = 4 – 17 years	Girls H = 75 – 180 cm A = 4 – 17 years
FVC	[l]	$H^{2.7}$	$0.95 * H^{2.7}$
IVC	[l]	$H^{2.7}$	$0.95 * H^{2.7}$
FEV1	[l]	$0.84 * H^{2.7}$	$0.81 * H^{2.7}$
MVV	[l]	$0.84 * H^{2.7} * 37.5$	$0.81 * H^{2.7} * 37.5$
FEV1%VC	[%]	84	84
PEF	[l/s]	$8.2 * H - 6.8$	$6.7 * H - 5.3$
FEF50%FVC	[l/s]	$5.6 * H - 4.4$	$4.6 * H - 3.3$
MEF50%FVC	[l/s]	$5.6 * H - 4.4$	$4.6 * H - 3.3$

Zapletal ⁸⁾		Boys H = 115 – 180 cm A = 6 – 17 years	Girls H = 115 – 180 cm A = 6 – 17 years
logVC	[ml]	-2.5768 + 2.7799 log(H)	-2.2970 + 2.6361 log(H)
logIRC	[ml]	-2.79590 + 2.73794 log(H)	-2.69813 + 2.67126 log(H)
logERC	[ml]	-3.81064 + 3.12550 log(H)	-2.74262 + 2.61668 log(H)
logVT	[ml]	-1.3956 + 1.8643 log(H)	-1.3956 + 1.843 log(H)
logFVC	[ml]	-2.9239 + 2.9360 log(H)	-2.7040 + 2.8181 log(H)
logFEV1	[ml]	-2.86521 + 2.87294 log(H)	-2.60565 + 2.74136 log(H)
FEV1%VC	[%]	90.6043 - 0.04104 * H	90.6043 - 0.0410 * H
logPEF	[l/s]	-4.37221 + 2.34275 log(H)	-4.37221 + 2.34275 log(H)
logMEF75%FVC	[l/s]	-4.01648 + 2.1541 log(H)	-4.01648 + 2.15414 log(H)
logMEF50%FVC	[l/s]	-4.21684 + 2.17719 log(H)	-4.21684 + 2.17719 log(H)
logMEF25%FVC	[l/s]	-4.58082 + 2.21169 log(H)	-4.58082 + 2.21169 log(H)
MVV	[ml]	-1.9178 + 3.0388 log(H)	-1.9178 + 3.0388 log(H)

GLI – Global Lung Function Initiative ⁹⁾	Men and women Age = 3 – 95 years (FEF25–75% and MEF25%FVC: 3 to 90 years)
---	---

The predicted values are calculated for:

- FVC
- FEV1
- FEV1/FVC
- FEF25–75%
- FEF75%FVC
- FEV075 (only children 3 to 7 years, Caucasian)
- FEV075/FVC (only children 3 to 7 years, Caucasian)

The predicted values are calculated depending on age, gender, height and ethnicity (African American, North-East Asian, South-East Asian, Caucasian and other/mixed).

The predicted values are calculated using this formula:

$$M = \exp(a_0 + a_1 \cdot \ln(\text{Height}) + a_2 \cdot \ln(\text{Age}) + a_3 \cdot \text{AfrAm} + a_4 \cdot \text{NEAsia} + a_5 \cdot \text{SEAsia} + a_6 \cdot \text{Other} + \text{Mspline})$$

For the coefficients a1, a2, a3 etc. there are defined value tables for each measured value of GLI from which the corresponding values are inserted into the calculation formula. The calculation formula remains the same for all the measured values mentioned above. Mspline is an age and measured value-dependent coefficient which is also taken from a value table defined by GLI.

The value PEF is not calculated when using the predicted value GLI.

Therefore, no predicted value curve can be mapped in the coordinate system in custo diagnostic.

Kainu (Finland)¹⁰⁾**Men and women**

Age = 18 – 83.99 years, Ethnicity: none

The predicted values are calculated for:

FEV1, FVC, FEV1/FVC, MEF75, MEF50, MEF25, MMEF (FEF25–75), PEF, FEV6, FEV1/FEV6

The predicted values are calculated depending on gender, height and age.

The predicted values are calculated using these formulas:

$$M = \exp(a0 + a1 * \ln(\text{height}) + a2 * \ln(\text{age}) + \text{Mspline})$$

$$S = \exp(b0 + b1 * \ln(\text{Age}) + \text{Sspline})$$

$$\text{LLN} = M - 1.645 * S$$

For the coefficients a0, a1, a2, b0 and b1, there are defined values for each measured value of Kainu,

which are inserted into the calculation formula to calculate the respective predicted value.

The calculation formula remains the same for all the measured values mentioned above.

Mspline and Sspline are age and measured value-dependent coefficients

which are also taken from a value table defined by Kainu.

Siriraj, Thailand¹¹⁾**Men**

H = 155 – 185 cm | A = 18 – 80 years

FVC	[l]	$-2.601 + 0.122 * A - 0.00046 * A^2 + 0.00023 * H^2 - 0.00061 * A * H$
FEV1	[l]	$-7.914 + 0.123 * A + 0.067 * H - 0.00034 * A^2 - 0.0007 * A * H$
FEF25–75%	[l/s]	$-19.049 + 0.201 * A + 0.207 * H - 0.00042 * A^2 - 0.00039 * H^2 - 0.0012 * A * H$
PEF	[l/s]	$-16.895 + 0.307 * A + 0.141 * H - 0.0018 * A^2 - 0.001 * A * H$
FEV1/FVC	---	$19.362 + 0.49 * A + 0.829 * H - 0.0023 * H^2 - 0.0041 * A * H$

Women

H = 155 – 185 cm | A = 18 – 80 years

FVC	[l]	$-5.914 + 0.088 * A + 0.056 * H - 0.0003 * A^2 - 0.0005 * A * H$
FEV1	[l]	$-10.6 + 0.085 * A + 0.12 * H - 0.00019 * A^2 - 0.00022 * H^2 - 0.00056 * A * H$
FEF25–75%	[l/s]	$-21.528 + 0.11 * A + 0.272 * H - 0.00017 * A^2 - 0.0007 * H^2 - 0.00082 * A * H$
PEF	[l/s]	$-31.355 + 0.162 * A + 0.391 * H - 0.00084 * A^2 - 0.00099 * H^2 - 0.00072 * A * H$
FEV1/FVC	---	$83.126 + 0.243 * A + 0.084 * H + 0.002 * A^2 - 0.0036 * A * H$

Danish Reference Values¹²⁾**Men**

H = 155 – 200 cm | A = 20 – 90 years

FEV1	[l]	$-2.87615 - 0.00026 * A^2 + 0.04201 * H$
FVC	[l]	$-5.17591 - 0.00026 * A^2 + 0.06015 * H$
FEV1/FVC	---	$105.77443 - 0.00126 * A^2 - 0.12261 * H$

Women

H = 150 – 195 cm | A = 20 – 90 years

FEV1	[l]	$-1.35015 - 0.00024 * A^2 + 0.02923 * H$
FVC	[l]	$-2.80132 - 0.00023 * A^2 + 0.04203 * H$
FEV1/FVC	---	$105.57449 - 0.00165 * A^2 - 0.12431 * H$

References for predicted values

- 1) Austrian reference values – sources
 - SPIROMETRY script of the Austrian Society of Pneumology
 Compiled by members of the Working Group for Clinical Respiratory Physiology, Standardization and Assessment
- 2) Cherniak – sources
 - PH. H. Quanjer et al.: "Lung Volumes and Forced Ventilatory Flows" Report Working Party, Standardization of Lung Function Tests, European Community of Steel and Coal; The European Respiratory Journal 1993, Volume 6, Supplement 16, 5 - 40;
 - "Standardization of Lung Function Tests in Paediatrics" The European Respiratory Journal, Volume 2, Supplement 4, March 1989, ISBN 87-16-14801-0;
 - R.M. Cherniak, M.B. Rober: "Normal standards of ventilatory function..." Am. Rev. Respir. Dis. 1972;
 - R.M. Cherniak: "Ventilatory function in normal children" Canad. Med. Assoc. 1962;
- 3) Knudson – sources
 - PH. H. Quanjer et al.: "Lung Volumes and Forced Ventilatory Flows" Report Working Party, Standardization of Lung Function Tests, European Community of Steel and Coal; The European Respiratory Journal 1993, Volume 6, Supplement 16, 5 - 40;
 - "Standardization of Lung Function Tests in Paediatrics" The European Respiratory Journal, Volume 2, Supplement 4, March 1989, ISBN 87-16-14801-0;
 - R. J. Knudson, M.D. Lebowitz, R.C. Slatin: "Normal standards variability, and effects of age" AM. Rev. Respir. Dis. 1983;
- 4) Ulmer – sources
 - PH. H. Quanjer et al.: "Lung Volumes and Forced Ventilatory Flows" Report Working Party, Standardization of Lung Function Tests, European Community of Steel and Coal; The European Respiratory Journal 1993, Volume 6, Supplement 16, 5 - 40;
 - W.T. Ulmer et al.: "Lungenfunktion: Physiologie und Pathophysiologie" Methodology, 2nd Edition, Thieme Stuttgart, 1983;
 - W.T. Ulmer, G. Reichel, D. Nolte, M.S. Islam: "Die Lungenfunktion" 3rd Edition, Thieme-Verlag Stuttgart, 1983;
- 5) Baur – sources
 - X. Baur: "Lungenfunktionsprüfung und Allergiediagnostik" Dustri-Verlag, Dr. Karl Feistle, 1998, ISBN 3-87185-270-8;
- 6) ECCS – sources
 - PH. H. Quanjer et al.: "Lung Volumes and Forced Ventilatory Flows" Report Working Party, Standardization of Lung Function Tests, European Community of Steel and Coal; The European Respiratory Journal 1993, Volume 6, Supplement 16, 5 - 40;
 - W. Schmidt: "Angewandte Lungenfunktions-Prüfung" 3rd Edition, 1987, Dustri-Verlag, ISBN 3-87185-130-2;
 - G. J. Tammeling, PH. H. Quanjer: "Physiologie der Atmung" Thomae;
 - H. Löllgen: "Kardiopulmonale Funktionsdiagnostik" CIBA-GEIGY GmbH;
- 7) Schindl – sources
 - R. Schindl, K. Aigner: "Atemfunktionsscreening und Sollwertebezug bei Kindern und Jugendlichen";
- 8) Zapletal – sources
 - A. Zapletal., M. Samánek, T. Paul: "Lung Function in Children and Adolescents" Progress in Respiration Research, Volume 22, Kager-Verlag, ISBN 3-8055-4495-2;
 - PH. H. Quanjer et al.: "Lung Volumes and Forced Ventilatory Flows" Report Working Party, Standardization of Lung Function Tests, European Community of Steel and Coal; The European Respiratory Journal 1993, Volume 6, Supplement 16, 5 - 40;
 - "Standardization of Lung Function Tests in Paediatrics" The European Respiratory Journal, Volume 2, Supplement 4, March 1989, ISBN 87-16-14801-0;
- 9) www.lungfunction.org
- 10) Reference values of spirometry for Finnish adults
 A. Kainu, K.L. Timonen, J. Toikka, B. Quaiser, J. Pitkäniemi, J. T. Kotaniemi, A. Lindqvist, E. Vanninen, E. Lämsinies and A. R. A. Sovijärvi
- 11) Siriraj, Thailand
 Dejsomritrutal W., Nana A., Maranetra N., et al. Reference spirometric values for healthy lifetime nonsmokers in Thailand. J Med Assoc Thai 2000; 83:457-466.
- 12) Danish Reference Values:
 Løkke A, Marott JL, Mortensen J, Nordestgaard BG, Dahl M and Lange P. New Danish reference values for spirometry. Clin Respir J DOI:10.1111/j.1752-699X.2012.00297.x.

3.8.3 Keyboard navigation and short-cuts in custo diagnostic

Use the quick links in the main navigation, the keyboard navigation and the keyboard short-cuts to enable fast and convenient working.

Quick links in the main navigation

User	custo med GmbH	1	
Patient		2	2
Examination		3	3

LEFT-CLICK

- 1 Change user's password
- 2 Call last patient
- 3 Examination main menu

RIGHT-CLICK

- 1 Evaluation search
- 2 Call last patient
- 3 Most recently opened evaluation

User	custo med GmbH	4	
Patient	Musterfrau Martina	5	4 10.10.1978 (39 Y.)
Examination	Spirometry	6	5

LEFT-CLICK

- 4 Change user's password
- 5 Patient master data
- 6 Menu for the current examination

RIGHT-CLICK
















- 4 All evaluations for the patient
- 5 Most recently opened evaluation for this examination

Keyboard navigation















When you press the Alt key, the initial letter of all the buttons on a screen page is underlined. Pressing an initial letter again triggers the corresponding button.

User	custo med GmbH	?	_	X
Patient				
Examination	Spirometry			
Holter				
ABPM				
Resting ECG				
Stress ECG				
Cardiopulmonary Exercise Testing				
Spirometry				
Telemetry				
Cardiac Rehab				
Mobile Cardiac Telemetry				
Blood glucose				
Ecg Streaming				
Task Manager				

Generally valid keyboard short-cuts

	End, cancel
	Confirm, continue
	Program information
	User main menu
	Patient main menu
	Examination main menu
	Patient master data for the selected patient
	All examinations belonging to the selected patient
	List of the most recently opened evaluations (same as clicking on the arrow button at top right)
	List of the most recently opened evaluations
	Evaluation search
	Waiting room list
	Device list
	Switch to Metasoft
	Create system report, service e-mail

Generally valid keyboard short-cuts in an open evaluation

		Unconfirmed report input dialogue
		Medication input dialogue
		Open comparison ¹⁾
		Open trend ¹⁾
		Open Print dialogue ¹⁾
		Open settings ¹⁾
		Open Options menu ¹⁾

¹⁾ Keyboard short-cuts will only work if the corresponding button is available on the screen page.



Manufacturer's contact details:

custo med GmbH
Maria-Merian-Str. 6
85521 Ottobrunn
Germany

Phone: +49 (0) 89 710 98 - 00
Fax: +49 (0) 89 710 98 - 10
E-mail: info@customed.de
Internet: www.customed.de

