



# 12-Channel PC ECG

custo cardio 300 | Resting and Stress ECG



## ECG device custo cardio 300 (USB) / BT-A

- Mobile use on all PCs in the network (USB, Bluetooth)
- Exchangeable patient cable: banana plugs, clips, D-Sub adapter for suction system
- ECG quality control directly on the device
- Start/stop button for manual ECG recordings
- Online analysis and display of heart rate related and ventricular events during recording
- Easy and cost-effective extension to CPET

*Device variants and type designations*


- USB version: custo cardio 300
- Bluetooth with Li-Ion battery: custo cardio 300 BT-A

## custo diagnostic (software)

- Comparison of all patients' examinations
- Comfortable reports with predefined text modules
- Clearly arranged A4 paper printouts
- Paperless medical practice: Mailing of evaluations and reports as email, fax, PDF and as evaluation file (.cst)
- Versatile connection possibilities into existing IT structures (HL7, FDA-XML, SQL, GDT, SCP)



Your authorised custo med dealer



## Technical data custo cardio 300 (USB) / BT-A

|                               |   |
|-------------------------------|---|
| ECG channels                  | 12  |
| Sampling frequency            | 1000, 2000, 4000, 8000, 16000 (only USB), 32000 (only USB) Hz   |
| Deviation                     | < 1,5%  |
| A/D converter                 | 24 Bit  |
| Input impedance               | > 50 MΩ   |
| Amplitude quantification      | 1,525 μV/bit  |
| CMRR                          | > 93 dB   |
| Defibrillation protection     | Electrical strength 5000 V (Recovery time < 5 s)  |
| Power supply                  | custo cardio 300 (USB): USB cable (standard connection to PC)<br>custo cardio 300 BT-A: Lithium-ion rechargeable battery  |
| Max. power consumption        | 1,5 Watt  |
| Operating time Li-ion battery | Bluetooth mode and 4 kHz max. 10 hours  |
| Charging time Li-ion battery  | max. 3 hours in pure charging mode  |
| IT connection                 | USB (cable length 3000 mm), Bluetooth   |
| Bluetooth frequency           | 2.4 GHz ISM frequency band  |
| Bluetooth range               | generally 10 m, depending on ambient conditions   |
| Patient connection            | patient cable 10-wire with banana plugs or clips,<br>D-Sub adapter  |
| Dimensions                    | 118 * 78 * 23 mm (L * B * H)  |
| Weight                        | approx. 130 g   |
| Patient leads                 | approx. 1050 mm (extremities), ca. 700 mm (chest wall)  |
| Classification                | Protection class II (USB version), device with<br>internal power supply (BT-A, with Li-ion battery)<br>class IIa, type CF |
| Underlying standards          | DIN EN 60601-1, DIN IEC 60601-2-25  |

Components for  
ECG recording:  
custo cardio 300  
Carrying case incl. belt  
Patient cable  
Disposable electrodes

Connection to the PC:  
USB cable A-A, 3 m  
USB universal power  
supply unit  
Bluetooth 4.0 USB Adapter  
USB extension cable A-A,  
1,8 m

Optional:  
D-Sub Adapter  
for suction system



True Wave®  
Filter  
Technology



### Resting ECG software functions

- Single beat analysis, cumulative complexes, measurement table, interpretation
- Recording of rhythm strips

### Stress ECG software functions

- Online arrhythmia detection with acoustic and visual warning function
- Predefined default profiles with the possibility for customised configuration
- Control of treadmills, sitting and supine ergometers
- Automatic recording of blood pressure, SPO2 and lactate measured values
- Summary of all examination-specific values in a report