

## **Manufacturer Information “CodeMeter”**

The licensing of custo diagnostic 4.x is managed by the software protection system “CodeMeter” offered by WIBU-SYSTEMS AG.

WIBU-SYSTEMS AG has reported a security weakness  
<https://www.wibu.com/de/support/security-advisories.html>

Generally, the operation of custo diagnostic takes place in closed medical networks, which cannot be accessed by external parties. In case your router, firewall, and the computers are protected sufficiently, the risk can be considered moderate.

WIBU-SYSTEMS AG made available an update for remedial action (CodeMeter User Runtime) under the following link <https://www.wibu.com/support/user/user-software.html>

In case of an update, all computers, on which custo diagnostic 4.x is installed (including PC with CodeMeter dongle), must be updated. This can be done by the user, or contact your regional custo med distributor. Generally, the operation can be done remotely.

10/2020

custo med GmbH  
85521 Ottobrunn