



Holter Software

Functional Overview of the Software Packages

		Software packages and recorders		
		Holter software light custo flash 501/L recorder	Holter software standard all recorder types ¹⁾	Holter software professional all recorder types ¹⁾
Software functions	Comparision of evaluations	✓	✓	✓
	Analysis (editing beat classes)	✓	✓	✓
	Trend/ECG	✓	✓	✓
	Examples for every event	✓	✓	✓
	Trend-Overview	✗	✓	✓
	Total ECG	✗	✓	✓
	Export as PDF file	✗	✓	✓
	Reduction of the data volume	✗	✓	✓
	Combination evaluation with ABPM ²⁾	custo screen 300/400 ³⁾ ✓	custo screen 300/400 ³⁾ ✓	custo screen 300/400 ³⁾ ✓
	multiday evaluation ²⁾ up to 7 or 3 days	✗	custo flash 500/510, custo watch ✓	custo flash 500/510, custo watch ✓
	Pacemaker ²⁾	✗	✗	custo flash 510 ✓
	Event system	✗	✗	✓
	ST Measurement	✗	✗	✓
	RR-Variability with scatter diagram	✗	✗	✓
	Heart Rate Turbulence	✗	✗	✓
ANS Diagnostics	✗	✗	custo watch ✓	

- 1) With custo holter software standard and professional all types of custo med holter recorders can be used, except for custo flash 501 and custo flash 501/L
 2) Software functions "combination evaluation", "multiday evaluation" and "pacemaker" are only possible with the recorder types mentioned in the table.
 3) custo screen 300 is an ABPM recorder. If it is attached to the patient simultaneously with a holter recorder the two evaluations will be displayed as Holter-ABPM evaluation based on their time-shared recording.

Your authorised custo med dealer