



custo screen protect hygiene set

Hygiene and comfort when wearing the blood pressure cuff



custo screen protect hygiene set:

- six custo screen protect fleece pads (with washing bag) for custo med blood pressure cuffs, sizes standard, X-Large and XX-Large

Utilisation:

- For each examination a clean custo screen protect fleece pad is applied underneath the blood pressure cuff
- The fleecy side is worn on the skin, the protective layer outwards
- The skidproof protection layer improves the fitting of the blood pressure cuff
- The fleecy inside makes for more wearing comfort
- custo screen protect fleece pads can be washed in the washing machine at 30 degrees – minimum 100 washing cycles



Your authorised custo med dealer

