



# custo sensitive

## the skin-friendly ECG electrode



- Suitable for holter, resting ECG and stress ECG
- Excellent ECG signal in combination with custo med ECG devices
- Ideal for stress ECG with a treadmill – excellent ECG quality, even at high speeds above 14 km/h
- Reduced motion artefacts due to central connection and solid gel technology
- Skin-friendly material, permeable to air and moisture
- Latex-free, PVC-free

Electrode size	approx. 50 * 35 mm (L * W)
Skin contact surface	approx. 43 * 35 mm (L * W)
Height without connection	1 mm
<b>Material</b>	
Sensor	Silver/silver-chloride (Ag/AgCl)
Gel type	Solid gel
Adhesive	Acrylic
Connecting part	Stainless steel snap fastener
Carrier foam	PE
Cover layer	APET
Label	PET/PP
<b>Shelf life</b>	
In opened bag	1 month
In unopened bag	24 months from the date of manufacture
Recommended storage temperature	15 – 30 °C
<b>Packaging</b>	
Bag	of 30 pieces
Box	of 40 bags (1200 pieces)
Smallest shipping unit	1 bag



Your authorised custo med dealer

