




# Medical ergometer

## custo ergo control 3000 e

Brake system	Eddy current brake rpm-independent Output range 6 – 1000 Watts adjustable in 1 Watt increments
Display	LC display for showing all important parameters LED display for rotational speed of proband display unit can be rotated by 180°
Adjustment options	Seat height from 75 – 110 cm, Infinitely adjustable grip (rotation: 360°) Height of the grip infinitely adjustable (60 – 120 cm) Adjustment screws to compensate for uneven floors
Dimensions	approx. 90 * 46 * (90 – 135) cm (L * W * H)
Weight	approx. 64 kg
max. patient weight	160 kg, 200 kg with anti-tilt mechanism (optional)
Power supply	100 – 240 V / 50 – 60 Hz
Interfaces	RS-232 for PC ECGs, USB
Protection class	II/B or II/BF (with blood pressure r
Classification	Ila in accordance with 93/42/ EEC
Applied standards	DIN EN 60601-1, DIN VDE 0750-2

- Ergonomic design
- Infinitely adjustable saddle and handlebar  
provide performance-enhancing seating position
- Highest safety due to stabile, mechanical construction
- High-quality material that can be easily disinfected meets  
highest hygienic standards
- Easy, intuitive operation
- Can be upgraded with automatic blood pressure  
monitoring module, SPO2 measurement, and powered  
seat height adjustment

Your authorised custo med dealer



 **custo-med**  
EXCELLENCE IN DIAGNOSTICS