


*Innovations in electrodes –
Developed by the leading ECG manufacturer*

Discover the secret of the perfect ECG!




NEW: *custo wing* – the new patch electrode for the perfect ECG!


No matter if Holter, resting or stress ECG - with the novel patch electrode *custo wing* the ECG will be perfect. Due to a hybrid gel, it is excellently suited for all areas of application. The decentralized snap fastener connection is fixed by an additional adhesive zone, so that artefacts caused by cable movements can be ruled out.



Our benefits:

- Hybrid gel for Holter, resting or stress ECG
- ★ Outstanding signal quality due to fully covered adhesive surface and decentralized connection
- ★ Artefact reduction due to adhesive, decentralized electrode connection as strain relief
- Large gel surface for reducing contact resistance
- Open-pored and breathable carrier material
- Very good initial and long-term adhesion (even with perspiration)
- Biocompatibility
- ★ Optimum storage due to resealable packaging
- ★ Fast and decent removal due to pre-punched carrier material





Suction electrode - ECG belt - Patch electrode *Innovations in electrodes from the ECG expert!*

2009: Suction electrode



With the exchangeable and reusable suction electrodes for custo cardio 200/400, custo med has introduced a completely new kind of electrode to the market and received a clinical certificate of hygiene for it.

2014: custo belt

This first medical ECG belt makes it possible to get a perfect 3-channel ECG without disturbing cables. It is washable, reusable, and highly comfortable for the patient to wear.

NEW: custo wing

The custo wing electrode now combines the benefits of the suction electrode and the custo belt. In combination with the new Holter ECG recorder guard holter, it provides an optimum, autonomous 3-channel Holter ECG. Thanks to the hybrid gel it is also excellently suited for stress ECG.



Innovative system solutions for your diagnostics


Since 1982, our brand has been representing innovative system solutions in cardiopulmonary diagnostics. In cooperation with committed and competent sales and service partners, we offer professional solutions taking account of your individual organisational circumstances and IT structures. No matter if your installation is based on a stationary or mobile system, a single-user workstation or multi-site network solution, in a doctor's office or in a hospital – your modular custo med system will be a profitable investment for the future.

Our modular platform custo diagnostic is the intelligent link between clinical engineering and administrative software.



Visit us
on YouTube!

www.customed.de



custo med – your local dealer

Contact:
custo med GmbH
Maria-Merian-Str. 6
85521 Ottobrunn
Germany

Phone: +49 89 71098 00
Fax: +49 89 71098 10
Web: www.customed.de
E-mail: info@customed.de