



2in1Suction System
and ECG



custo cardio 400/400 accu

12-channel resting and stress ECG with integrated suction system

custo cardio 400/400 accu sets new standards in functionality, usability and hygiene. Thanks to the exchangeable, reusable suction electrodes a hygienic and sustainable application is guaranteed. The version with rechargeable battery is completely self-sufficient and can be wirelessly integrated into your WiFi network, making it flexible to use.



custo cardio 400/400 accu*:


- · Quick and uncomplicated cleaning
- · Perfect suction capacity in six steps
- · Stable and fast data transfer
- No room dependency, can thus be applied with full flexibility and mobility*
- · More than 100 resting ECGs with one battery charge*

Different electrodes/suction lines:

- · Gentle application due to adapted pressure
- · Electrodes and suction lines can be easily exchanged
- Cleaning of suction lines by automatic blowing after each examination







Resting ECG Highlights in the custo diagnostic software

The complete ECG at a glance!

ECG overview of all channels over the total recording period. The area to be analysed is automatically adjustable or freely selectable.

All important values are visualised compactly.

Comprehensive measured value tables support the automatic measurement and interpretation.


Useful additional options.

Options such as HRV analysis and sport ECG enhance the resting ECG, as does the option of an evaluation comparison which compares two evaluations in detail.

Time-saving and precise reports.

With the help of the individually created text modules, the proposal is quickly created. The printout summarizes all results and the proposal clearly on one page.

Versatile ECG interpretation program.



NEW: Sport ECG as an option

Stress ECG Highlights in the custo diagnostic software

Easy control of the examination.


Uncomplicated control of treadmills, cycle ergometers and recumbent ergometers. The module provides predefined standard profiles, with the possibility to configure them freely. Even during a running stress test, manual intervention in the load profile is possible.

The patient is continuously monitored.

Beside a detailed documentation, there are a variety of notification options with acoustic or visual signal. In addition to the ECG lead, the measured values of blood pressure, SpO2 and lactate can be recorded automatically. Of course, the ST segments are measured continuously and displayed graphically.

Time-saving and precise reports.

With the help of the individually created text modules, the proposal is quickly created. The printout summarizes all results and the proposal clearly on one page.



Clear total view of all important values

Since 1982, our brand has been representing innovative system solutions in cardiopulmonary diagnostics. In cooperation with committed and competent sales and service partners, we offer professional solutions taking account of your individual organisational circumstances and IT structures. No matter if your installation is based on a stationary or mobile system, a single-user workstation or multi-site network solution, in a doctor's

office or in a hospital – your modular custo med system will be a profitable investment for the future.

Our modular platform custo diagnostic is the intelligent link between clinical engineering and administrative software.

Visit us on YouTube!

www.customed.de

ingle-usc. Resting Stress **ECG**

Central

Blood

Pressure

Cardiac Rehab

Spirometry

Hospital Information System

PACS

custo med - your local dealer

Contact: custo med GmbH Maria-Merian-Str. 6 85521 Ottobrunn Germany

Phone: +49 89 71098 00 +49 89 71098 10 Web: www.customed.de E-mail: info@customed.de

