


custo screen check kit



With this kit you are able to check the accuracy of all custo screen recorders according to norm EN 1060-1, 1060-2, 1060-3. The screen check software is suitable for testing all the long-term blood pressure measurement devices screen 300/400 and also the older screen 100/200 devices. The screen check software runs under Windows OS (Win XP, Win 7, Win 8, Win 10)

The kit contains the following items:

- Solid suitcase with inlay
- One year license "custo screen check"
- Code meter USB key with the custo screen check software.
- Instruction manual and video on USB stick
- Infrared custo com IR communication interface
- Dead volume 500ml
- Hand pump and connection tubes
- Calibrated pressure measurement device

For more information contact service@customed.de